

Product datasheet for **SC214482**

ARL6IP6 (NM_152522) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ARL6IP6
Synonyms:	AIP-6; AIP6; PFAAP1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_152522
Insert Size:	1437 bp



Insert Sequence:

>SC214482 3' UTR clone of NM_152522

The sequence shown below is from the reference sequence of NM_152522. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site **Blue**=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

GAATGGCATCGTAGCTGCTCTTACTGTAGCATGGTGCCTCATG**AA**ACCCCACTGGAGCGATATTGTTG
GCAAACTTAATCATGATTGTTTTGTAATAACAAGAAGGAGCATCACTGTCTACTCAGAAGACTGAGAAA
CCTGCTGTTTATTATGTAGTTCAGATATTTATCACAATCATCTCATTATGGAAGACCTTTTAAAGCATT
GTTTTAGAAATGCTGAGTATTAAGATACAGATTAATTGGGAATATCTGAGTATTAAGTACTTTCTTCAGA
GTATAAGATGTTTACCTCATCTTTTACTTTTGTGTGTAGTTCTTTCAAGTTGTAGGAAACATTTAA
TGGAATTTAACTCAAACTTGAATACAGGACAATGCTTGCCTTTTCATGTATGTACTACATTTTTTGC
TAAGAAATACTGATATTTCTGTTTAGTTGAGCTAGAAATACTTCTTATTTATACATTTAGGAAAGCAAAT
AATGCCTACTACTCCGACTTTTATAGAAGCTACTTTTAAATCAGAATATTTAATTTTTGATATTCATATA
ATTAATAGAAGTTCATTTATATTTTTATGGGGCATAGTTCCTTATGTGTTTTTAAATGTATTTTCAT
ACTACATACTGAATTTGTATGTTTTTAAATTTGTACATCTAGACAACGTAAACATTTATTTTTAGCT
AGTGCAAACTAGTACCTGCCATTTTACTAATTTTGTCTTTAAAAAAGCAAAAAGCACATTGACCT
AAGTTGAAAATTAAGTAAGTTTATTTTTAACAAAAATGCCTGAAAAAGCTAAATTTTAAAGTCATTA
AGATATTGAGAAAAAATTTGAAATTTTACTATTTCTGTTTCCACAATCCAAATATTTATCTTGGTGTA
TATATTGTTACTTTAACAGAACTTGCAATTTTGTGTTTTAATAAATATAAACATGAAATTTTGTATGT
GAGAATGATTGAAGTAGTTTGTCTTAATCTCAAAATTTAGTTACCAAGTAGAAAAGGTATTTTGATA
CTAGATATTTAACTACATATAGTTAATATAATTTTATAATTTTGAATTAATGTGCTCCTGTAATT
GTGATATTTTGTGTTCTAAGTAATTAAGCCTTAATTTTCCCTATGTTACTAAAGACTTTTAAATGTTTA
GAAAGTTACCTCAGTTTGTAGAAAGATGGACTACTTACAAAGCTGTTTTCCCTGGCCATAGGAAAACCTTAC
AATAAGAAACCATTAAGTAGTCAAACTGCTTATTCAGGCTAAGTGGGAAATGAATTTGTGCAGTGTGTT
CCTTAATTATCTACTGAATGTTGTTACCGACTAAACAAGGTTTCTAAAAGTCTCGCATTTCTTTTACTAT
TCAAACTCTTTGGTTGGAAATGACTAATACTGCTCC

ACGCGTAAGCGGCCCGGCATCTAGATTCAAGAAAATGACCG

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_152522.4

Locus ID:

151188