

Product datasheet for **SC214509**

FBRSL1 (NM_001142641) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FBRSL1 (NM_001142641) Human 3' UTR Clone
Symbol:	FBRSL1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001142641
Insert Size:	1440 bp



[View online »](#)

Insert Sequence:	>SC214509 3'UTR clone of NM_001142641 The sequence shown below is from the reference sequence of NM_001142641. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CGAAATCCCCCGAGGTGGAGGCGCGGTAGCCCCGGGGCCGAGACGCCTCTCCGAGCGGAGCGCACCG CTGTCCGTCTCTCCATCAGTTCCCTAGAACTCAAGCACAGCTCCCGCCGATCCTGGGGCGGCGCCGCTC TCCACCCGAGCCTGCGGAGGCGGGGACTTGGGTGTGCGCTTTTCTGAGGGTGATCTTTTGTTCGGA GTTTGGGATTCGCTGTTGGGAAAAAAATCGCGTTTGTACCTTTTCCCCCACAGATGAGAAAGTGT TGTAATGGATTTGATTTTTCTTAATTTTATGCTCTTTAGACGTTTAAAGAACCAAAAAAGAACTTT TGCTTCTTCGTTTTCGTTGGGATTTGCAGTTTAAATGGCCTCAGATCCTTCTCCAGATCCCAGGGTTT CTTTGTCTTATTTATGGAGAAAAACCGTCACTTTGTCCAGCGCACTGTGAGGCCCCCACTCAGGCCAG CCCTGGCCCCCCTTGGTACTTGAACCGAAGTTACAGATTATATTAATAAATAATGTACAAAACCTTT GTTTTTTTGCCTTATTATGCAGAAGTGTAATTTCTTTTTTGGTTTGTTTTTTAAAGCAACGAAACG TGTAATAGTCTGTTGATATAAATATATGACGTTACTAAATCTTAATCTAGATAGACTTTATAAAAC CGTTTCCAGAAACCCCTCTGTTGTTTGTCTAACATGGTCACAAAAACCCAGGCGCTGCTGGCTGTAAA CATTATCAGAAGTTAATGGCAGCACTTTCTTCAACTATGCAAGCGCTTCCCGCGGCTCTGCTGGT TGGGGGAGGAAGTTCCGCCAGCGACCTCCAGCCCCGTTGGTCCGTGGAGGGTTGGTCCCCCAGCAG AGCCCCCTGGGCCAGCAGAACAGCACTCCTGGGCGGCTCCTGTGGGCTTCAGGACTGGCCGGCTCCAGCC GCAAGGGCGATGCTGTGGACGTTTGTGTGGATGGTGCAGTGTTCTTGAATCTGTAGACCAGATCGGTG GAGGTATAGAAGTGCAAGAGTCTAAAGGCTGATTTACACACACACACACACACACACACACACACA AAATCACTGCTGTGTGTTTTCTTCCGCTGAGGCATCAGCGTCACCTCGGTGTGGTGTGCTCGCTCCAG GCCCTTGGCCTGGCCTCCTCTCAGTTTGGGTTTTGCCCTCACCCACATGAGCCCCCATGCCCTGCC CCCTTGGCGCCTTTTGTGTTCCCGATGTCTGGGCTGTGTCCTGCCGGGATTCCAGCCTTCGGGGTC AGCCGTCGAGCCACTTGCTTTGCTCAGTTCTGTGTTTCAGCCACATTTCCAGGTGGACCCTGTTTGT TTAAATATTCTGTTCCCATGTCAATCAGTGACGATAAATACAGCCTTGATTTGGATGAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001142641.2</u>
Locus ID:	57666
MW:	54.7