

Product datasheet for **SC214514**

ZFP28 (NM_020828) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZFP28
Synonyms:	mkr5
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_020828
Insert Size:	1446 bp



Insert Sequence:

>SC214514 3'UTR clone of NM_020828

The sequence shown below is from the reference sequence of NM_020828. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TGGAAATCCATCCTCCCTCCCATCACCATAGCCTCGAGACGTCATTTCTGTTTGACTACTCCAGCAGTTT
AAAACCCCATCTCCCTGCCCTTTTGTCTTTTGTCCCTTATTAGTTAGTTCTTCACATAAGTGTA
ATGTAACCTATTCACTCCTCTGTGAAAACTTATAGTTTCTTTAAATTTGGTTAATGTGTGAGATGTGCTC
AGCACAGTGCCTGGTCCATAGTAAGTGCTCAGTAACTTAGCTGTTTTAAAACTTTGTATTTGAACAT
TGAAAAAGTTACAGTAGTCAGCTCTGATAAAAAATGATGCAGTAGGGTGAGGGTAGGAAAAAGCACATT
TTCTATCAGGAACAGAATTCTCCAGTAGTGGGTGAGGTTTTGCCTTTGTTGGTTTTAAAACTTGATTCT
ATAATGCCAAGTTAGTTTTGTGGCTTTCCATCTGACCCTATGTGAATGAAGGTGATGTGACCTTGTG
GTGAGAAAAATTAACCTTTACATTTGACTTGATTGTTTTAGAAAGTCTAGGGACCATGAATGAATAGGC
CAGCTGGGACAAATGAATTTAAAAATCAGAAAAATGCAAGATTTATATGCATGAAGTTAAACAACCTG
ACGTTACTCAAGAATTAGAAAACTTTGCAAGATTTGACTTGTTTTAAAAATCACATTTATAAGTGAACCG
TATTAACCTTTTAAGGAACCATTCATTGTGAGGTAACTGATCCAGAATAGGGGTGAGCAAACTATGA
CTCATGGCCACAGTCTCACTGACTGTTTTGTATGGTCCATGATCTAGAATTTAAAAAATTTTAAAGG
GTTGAAAAAAGTGAAAAAGAAATTTTTCAACATGAAAAATTATGAAATTTAAGTTTTGGTGTCTGTAAA
TAAAGTTTTGTTGGCTTCAGCCACACAGTGATTTACACCTTGATTGTGCTGCTTCCAGGCTGTAGCAGC
AGAGTTGAGCAGTTGTGACAGGAGACCATGTGGCCTGCAGAGCCCAATATCTACTGTCTGATCTTACA
CAGAAAAATGTGTGATCCCTGCTATGGAGCAGAGGTTATCAAACTAAAGCCCATGGACCAATCCTAT
CTGCTGCTTGTGTTTTGTAATAGAGTTTTATCAAACACAGCCATGCTTACTTTTGTAGCTATTGACTA
TGGCTGCTTTAGACAACCTGTGACAGACTATATGGCTCGAAAACTGCAAAATTTTCTATCCTTTAACA
GAAAGTTTGCCAACTCTGCTCTAGAGTAGAGAAAAATGTATAAAAGATTTTATTTTATGAGGGCAAT
ACAACTGTCACATCAGAAACAAGAAAAAATGATAAAGGAACTCATTTATCAATAGAGGTGAAAGGAA
ATTATTAACTATATTGAAAAATAAGATGTCAATAAAAGGAGAAATGATATTTTCTAGATGAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTGCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_020828.2](#)

Summary:

May be involved in transcriptional regulation. May have a role in embryonic development. [UniProtKB/Swiss-Prot Function]

Locus ID:

140612

MW:

55.7