

Product datasheet for **SC214528**

CCDC43 (NM_144609) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	CCDC43
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_144609
Insert Size:	1446 bp



Insert Sequence:

>SC214528 3'UTR clone of NM_144609

The sequence shown below is from the reference sequence of NM_144609. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGGACACAGAGAGGGGAGCGAAAGCGATACCTTGGCCTACTCTATTCTTTCCCCCATGAACTTATCA
AAGGGAGAAGTGGTTGTGTATTTAAGAGAAATGAGAGTGCATTGCAAGTAGTTTCCCAAGGCATTTCC
CCTATTGTTACATGGTCAAAAGGAAAAATCAAAATCACTTCTAATTACAAAAATGTGCTGTTTGGTGGG
GTGGGCAATCAGATTATTATAGTTGATGACTGTACCAAGATCTGGAATGGGGCCACCTTTATATTTT
AATACTCAAGTCCTGTGCTCCTTTTGGCCACTTTAAAAAGTCTTCCCTCACCTTGGATAAAATAGGATTC
AGAAAGCTAACAGTTTTTCACTTGAAGGAGCAGAAACATAATCAATCCCTCCCTGAGATAAAATGTGGG
CATGCACATTGCTTTGGAGTCAGTATTCAGCTGTAAGTGTAGTATTGTTTATGGCTGAAATAGAAAGC
TAATGAAGAGGCCCTACCTGGATCCAGATTTTATGTTAACATGGCAGCAAACTAGGAAAATGGGAAGAA
ATTACTTGGGGGTGTTCTAGATCTGATTTTTTTTTTTTTTTTGAAGTTTCAGCCTGTGGGAGGAAAA
GAGGTATCTTTTAAAAATAGAGATAGTTACATATTTTCTCATCTAGGCTGTGTCACTATTGGGTGCAAG
GATATCATTGTTGGCTAGATTCATTTTGTATAATCATGTATCCTCTTGTGTGCTGGTAGAGATTTAAT
CCTGATTTTTCCATAAAACATGAGTATTAAGAAATAATCCTGGTTTGGAGAACTGGAGAAATCACCC
TTTTAAGGAAGAAACACTGGAATTTCTGCTAACACCAAGATATTTAAGAGTGTACATAGTAGGTGCTC
AACAAATTTATTGAATGAATGAGTGAATGGAAAACTGGGAGAGTCAAAAGTGAGCAGAAGCTCTCCAT
TTCTACTTCTGTCAAAACCATTAATTTGTAATAAGGCCCTTCTCCACTTGACTTCAGGCAGCAGA
TTGTCTAGAAGCCTAAGGACAGCAATTTCTCTGACAAGACAAAGTAGATATTTTATACCAGGGGTGGC
AAACTACTGCCACGGGCCGAATTTGGCCAGTCTGTTTTGTATGGTGCAAACTAAAAATGATTTTA
CATTTTTAAAGAGTTATAAAAGAAAAAATATGTGGTCTGTGAAATCTAAAATATTTACTACCTGGCCT
GTTGGAGGAAAGGTTTGCCAATCTCTGTTTATACCATTAAGTATGAGATTAACAAAACTTTTACCTT
TGTGCAGAAGGTTAAAAAATACATGGTTAAGGAAAAGGAGACATGTTACCTCTTCATACACTCC
TATAACTGTGGCATTGCAAAAAATAAAAAATAACCACCTTTAAAAAATAAATCTTATTTAAATTGC
ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_144609.3](#)

Locus ID:

124808

MW:

55.8