

Product datasheet for **SC214537**

ZNF323 (ZSCAN31) (NM_145909) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF323
Synonyms:	ZNF20-Lp; ZNF310P; ZNF323
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_145909
Insert Size:	1459 bp



Insert Sequence:

>SC214537 3'UTR clone of NM_145909

The sequence shown below is from the reference sequence of NM_145909. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CAGAAAATCCACACTGGAGAGAGACCATAGGGTGATGAGTTTGGGAAAGCCTTCAGTCATCATTGCAA
CCTTATTAGGCATTTTGAATCCATACTGTTCCAGCAGAACTGGACTAATTCTGTGGACCCTCAGACCC
CTAAATGGGGCCGCTCTTGGAAAGTAAGATTTTCAAGTGGCTTATTCTGTTCTTATAAAATTTATTTTGT
ATTGCATTGAACCACTTGTCTACAAACCAGGTGTTGGTGTCTCCTGTGTGGATGACTGTGATTGTAGGG
GGTTGCTGGGTGGCCTTCTTATTCCTAGCCTTTGGTCCATTCTTATGATTTCTTTAAAGGCACCAA
ATCTTATCTGGAAATAAATTGTTATGGAGGGGCTGTTTTTGAGCAATACTGCATATTGAACAATTCTAA
AAGACCAGTATCTATGTATACTTTAAAAAAGATTGCTAAGATCAAAGGATTCAGAGGTGCATC
TTTTCTGTAAGAGTGACATTTTAAACAAATGTAAGGTTTTGGGAATAAGAATCGTCCACATCCCAC
AGGAAATCATTGCTATGCTCTCCCTAGGAATGTGATAGATAGCTTCTATAATTATGCAGCATTATCCA
TAGTATGAAACCTGTAGTTTTAATGGTTTGAATCTATGTCTCATTTTCATATCTAACTCATGTTTTTG
ATACGAATCAATATATATGTTGATATCATGAGTTAATATGAGTTTGACATGAATCATCAACATAGTATA
TACATATCAACATAACTATGTTGATGATTATCAAAATCTCAATGCAGGTCATTACTGTAAGTTTTTC
CATTTTTGTTACAAATAATAAGTAAGATGATACTTTAAATCTCTGTTCTTCAAGGTACCTCTTTAA
TTCATAGCCTTTTATGTCAGCATTATTCTTCTTTGGGGTGGTTCAGCAACTTTGAGTGCCATTGCAA
GCAGCATCATTTAAGTAGTTTTATAAGAAAAGAAGTCTGATATCTTGGCAGCTTGCTGCTTACTGGG
GAGACCAGCATGCTTAGACACTTATTGCCAGGGTTAAACATATGCTCTGGCACCAGTTCAAATCTTTC
CTCCCCATTTGTGAGCTGTGGACTTTAGGCAAGTTATCTAAATGCTGCTATGCTATAGTCTCCCCACC
TGTAACAGAGGATGATACCAATACCTTCTCATAGGGTACATAAGATTAAATAAATACAACTTGTA
CTTAAATGTAGCTCTCAATAAATATGCATTATGATTACATACGAAATAATAGCATTTAACAAATGTA
GGTAGAATTTTAAAGTATCATGTTAAATAGGATACGAAATGCCTCTGCTGACATTCTGTACTCTTT
TCTCCACAGTCAAAACCTGAATCAGAATATTTCAAAAGAGAGGAAATAAAGTTCAGCTTAAACTGTT
CTGCCCTAAA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG

```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_145909.3](#)

Summary:

This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]

Locus ID:

64288

MW:

56.6