

Product datasheet for **SC214540**

ZNF323 (ZSCAN31) (NM_001135215) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF323
Synonyms:	ZNF20-Lp; ZNF310P; ZNF323
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001135215
Insert Size:	1459 bp



Insert Sequence:

>SC214540 3'UTR clone of NM_001135215

The sequence shown below is from the reference sequence of NM_001135215. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CAGAAAATCCACACTGGAGAGAGACCATTAAGGGTGATGAGTTTGGGAAAGCCTTCAGTCATCATTGCAA
CCTTATTAGGCATTTTGAATCCATACTGTTCCAGCAGAACTGGACTAATTCTGTGGACCCTCAGACCC
CTAAATGGGGCCGCTCTTGGAAAGTAAGATTTTCAAGTGGCTTATTCTGTCTTATAAAATTTATTTTGT
ATTGCATTGAACCACTTGTCTACAAACCAGGTGTTGGTGTCTCCTGTGTGGATGACTGTGATTGTAGGG
GGTTGCTGGGTGGCCTTCTTATTCCTAGCCTTTGGTCCATTCTTATGATTTCTTTAAAGGCACCAA
ATCTTATCTGGAAATAAATTGTTATGGAGGGGCTGTTTTTGAGCAATACTGCATATTGAACAATTCTAA
AAGACCAGTATCTATGTATACTTTTAAAAAAGATTGCTAAGATCAAAGGATATTCAGAGGTGCATC
TTTTCTGTAAGAGTGACATTTTAAACAAATGTAAGGTTTTGGGAATAAGAATCGTCCACATTCAC
AGGAAATCATTGCTATGCTCTCCCTAGGAATGTGATAGATAGCTTCTATAATTATGCAGCATTATCCA
TAGTATGAAACCTGTAGTTTTAATGGTTTAGAAATCTATGTCTCATTTTCATATCTAACTCATGTTTTG
ATACGAATCAATATATATGTTGATATCATGAGTTAATATGAGTTTGACATGAATCATCAACATAGTATA
TACATATCAACATAACTATGTTGATGATTATCAAAATCTCAATGCAGGTCATTACTGTAAGTTTTTC
CATTTTTGTTACAAATAATAAGTAAGATGATACTTTAAATCTCTGTTCTTCAAGGTACCTCTTTAA
TTCATAGCCTTTTATGTCAGCATTATTCTTCTTTGGGGTTGGTCAGCAACTTTGAGTGCCATTGCAA
GCAGCATCATTTAAGTAGTTTTATAAGAAAAGAAGTCTGATATCTTGGCAGCTTGCTGCTTACTGGG
GAGACCAGCATGCTTAGACACTTATTGCCAGGGTTAAACATATGCTCTGGCACCAAGTTCAAATCTTTC
CTCCCCATTTGTGAGCTGTGGGACTTTAGGCAAGTTATCTAAATGTGCTATGCTATAGTCTCCCCACC
TGTAACAGAGGATGATACCAATACCTTCCTCATAGGGTACATAAGATTAAATAAATACAACTTGTA
CTTAAATGTAGCTCTCAATAAATATGCATTATGATTACATACGAAATAATAGCATTTAACAAATGTA
GGTAGAATTTAAGTATCATGTTAAATAGGATACGAAATGCCTCTGCTGACATTCTGTACTCTTT
TCTCCACAGTCAAAACCTGAATCAGAATATTTCAAAAGAGAGGAAATAAAGTTCAGCTTAAACTGTT
CTGCCCTAAA
ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_001135215.1](#)

Summary:

This gene encodes a protein containing multiple C2H2-type zinc finger motifs. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016]

Locus ID:

64288

MW: 56.6