

Product datasheet for **SC214571**

TBCID9 (NM_015130) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TBCID9
Synonyms:	GRAMD9; MDR1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_015130
Insert Size:	1443 bp



Insert Sequence:

>SC214571 3'UTR clone of NM_015130

The sequence shown below is from the reference sequence of NM_015130. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GACTATGAAATCTCGGCCATGTCCGGCTGACACGGGCGCCTTCCCGGGGAGTGGGAGGAGAGGGAGGG
GAGGGATTTTTTATGTTCTTGTGTGTGATTTTTTCTTTCTTTTAAATTAATATTTATTAGTA
CCTGGCTTGAAGCCTAGTGTTCATAATGTAATTCAATGAAACTGTTGGAGAAATATTTAAACACCT
CAATGTAGGTACATTACACTCTTGTGCGGGGAGGGGATTACCAGAATACAGTTTATTCGTGAATTC
TAAAAAACAAAAAGATGAATCTGTCAGTGATATGTGTATTATAACTTTAATCTTGCTGTTGAGCT
GTATACATGGTTTTAAAAATAGTACTGTTTAATGCTAAGTAAGGCAGCAGTCATTTGTGTATTGAGGCT
TTTTAAATAAAATTAGAGCTGAAGGAAATGAAAAGCCACAAATGCAAGACTGTTCTTAATGGAAGG
CATAGTCAGCGAGGGTAAATCCTATACCACTTTAGGAAGTATTAATAATTTTTAAGATTTGAAATAT
ATTTATAGAAAGTCTCTATTCAAAATCATATTCCACAGATGTTCCCTTCAAAGGGAAAAACATTTGGG
GTTCTAAACAGTTATGAAAGTAAGTGATTTTTACATGATTCAGAATAACACTTGTATTGACCAATTTA
GACAGATACCAGACCAATTTTGCAATTAAGAAATTGTTCTGATTATTTACGTCAACTCATTAGAATTCA
GTGAAAGTAACAGTCTTTTGTACAGAGAATCTGAAAGTAGCAGCAAAGACAGAGGGCTCATGACAGG
TTTTTGCTTTTGTCTTTGCTTTTGTGTTTTGAAAGAGTAAAGTACTGATGCTTCTGATACTGGATGTTA
GCTTCTTACTGCAAAAACATAAGTAAACAGTCAACTTTACCATTTCCTGATTCTCCATAGATTGAAGA
AATTTATACCACATATCGCATATGACCATCTTCCATCAAATCAATGTAGAGATAATGTAACTGAAAA
AAAATCTGCAAGATAATGTAAGTGAATGTTTTAAAAACAGAACTTGCACTTTATATAAAGAATAGTA
TGCTCTATTTCTGAATGGATGTGGAATGAAAGCTAGCGCACCTGCACTTTGAATTTCTGCTTCTTT
TTATTACTGTTATGATTTTGTCTTTTACAGATGTTGGACGATTTTTTCTTCTGATTGTTGAATTCATA
TCATGGTCTCATTTCTTTGCTTCTTTGGAATATTTCTTTCAACACATTCTTTATTTTATTATACATT
GTGCTCTTTTTTGTCTATTGCTGCTGTTGTTTTTATTCTATTACAGGATGATTTTTAACTGTCAA
ATGAAAGTAGTGTTAACCTCAAATAGGCTAAATGTGAACAAATAAAATACAGCAAATACTCAGA
ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_015130.3](#)

Summary:

May act as a GTPase-activating protein for Rab family protein(s).[UniProtKB/Swiss-Prot Function]

Locus ID:

23158

MW:

55.9