

Product datasheet for **SC214581**

USP9Y (NM_004654) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	USP9Y (NM_004654) Human 3' UTR Clone
Symbol:	USP9Y
Synonyms:	DFFRY; SPGFY2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004654
Insert Size:	1453 bp



[View online »](#)

Insert Sequence: >SC214581 3'UTR clone of NM_004654
 The sequence shown below is from the reference sequence of NM_004654. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GTATCCTCACCTCAGATGAAGGATCAGTGAAGCAATAATTAAGTCTTCTTTATGACTATGCACTA
AGGTCTTATAGTCCAACTTTCTCTGTGTCTGGCTAGTATTGAAACTAGATAAACTGCTCCAAACCAA
CATGGAGTAAAGAGCATATTTCACTGGTTTATTTGCAGTAATTTGCAATTTGTAGTGATAAGACACAT
GCAGGGTGAAGTACAGAGTTTGTAAACAAATGACTGGTCTAATCTGTAAATGAGAAAGGTATATAT
ACTATGTTAATGTCTGACTGTAAATTTCTAAGCAAGAACTTTTTTGTGAAACAAGTCAGATCTAC
ACAGTCACACAATTATTTTTGTGTGTTCACTACATTGTCAATTGATATTGCCTGCTTTGAGCAGTT
TGGTCAACTTACCAACTTCCCCCAAAAAAGGGAACATAAAAGAGCCCATCTTTGTAGTTTACACCAA
TAGTTTCTGTAAATCTTTCTCTGGATATATAAGGCTGGTGGTAACCTTTGAATTATATGGTTGAT
GTGGAAAATTGGCAGTGAACATTTCTAGATACTTTTCATTACCTTTTTATTCTGGTATATAGGCTAAC
CACTTTAAAGCTATTCTTATGCTGTAACAGTTAGCATGGCTTCACACTGTTTGTGTAGCCAAGAGGACA
GAATTACATGAATGACAGTGCCAGAGTGACAGCTGTATATTGCTCAGAGCTTTTATTTCTTATACCTA
GAATAATATAAAATGGGGGAAAAATGTGACAGACAAGCAGTTTTTCATTGCACACATATTTCTCACAT
TTAATGTTTGATAGTTCAATTTCTTGAGCTGAAGAATATCAAAGTATCGATAATACGAATTTTAAATA
TTAAACCAATACACAAAATTTCTATGTCAGAAATGTGGTGGAGCATAATAGATTGTATTTGGTGTGC
TTGCGATTTTTTTTTCCATAGAATTTATTAAGTGAAGTTCTAAACTTTGCTTCTCTGATCCCGGT
GAAGTGATACATCATAAGAATCCATAGTACTTTGAAGTACCATTGCACCAAGATGTCTGACTGAATTCAT
AGTCACACTTTTTATTTGAAAGAAAGAATTGTTGTAGTTTTTTTTTATTATTCTAAACTCTTGTGTGA
GATACAAGATTTAATTAAGATCTAAGCTCTGCTTATTTAATGTAATTCTAAGGTACCATTTTAGAAAA
AACATTTGTTTTAAGATTTCCAAGAACTGTGAGTTAATACTATATTTAAAGAGAATTGGTAAATTTT
GAATGTGTGAATATTTTGAACCTGTTTAAAAACCAATATACCTGCAATAGATACAGCCTATCCTA
TACTATTTAAATGTTGGCTGTTTTGTTTATAGAAATATTTTGTGTAATTCACAAATAAAATTTAAG
AAGA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_004654.4](#)

Summary: This gene is a member of the peptidase C19 family. It encodes a protein that is similar to ubiquitin-specific proteases, which cleave the ubiquitin moiety from ubiquitin-fused precursors and ubiquitinated proteins. [provided by RefSeq, Mar 2009]

Locus ID: 8287

MW: 57