

Product datasheet for **SC214598**

PSF3 (GINS3) (NM_001126130) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PSF3 (GINS3) (NM_001126130) Human 3' UTR Clone
Symbol:	PSF3
Synonyms:	PSF3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001126130
Insert Size:	1448 bp



[View online »](#)

Insert Sequence:	>SC214598 3'UTR clone of NM_001126130 The sequence shown below is from the reference sequence of NM_001126130. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AAGAGAAAATTCAGTATGGAAGACBAAAGCCGGAAGAACACAGAATGGCTCCTCACAGACGTATC CCTCCGTGTGCTTGATAGGAGCTGGTTGACCTTGACAGAACCAAGATCCTGTCCCATTTTCATGGCT TATTTCTGTGGCCATAGAGAATTATAGGGAACGGACATGCTGGAGGATGTGGGTGTCCCTGGCTCTG TGAGTCTTCCAGGACCGTCCACCCTGCTGACCCACAGCCAGGCCCTTTAACCAAGAACCCATGGCC AAGGAGAAATCAAAGTCCTTCTAAATAAGAATCACTGCCATATAATATACAGTAGAGTTGCAACT GAGATTCCTTGTGTCTGGGAGTTTGGACAGCTTCAGATGTACAGTTTCTAGCCACAAAGCACAGGTA CAAAGTGGGTATCGCTGTTCAAAAAATGCTCTCTTGATCTTATTTGCCTCATCTTCTCATGGTTGT ACAGAGGATAGCACCCACCATGCCAGCCTGACTTGGAGATATCTCTGCTGCCTGCCTGCAGGGAGTT ACCCACAGTTTCAAAAACAGTCGCCAGATAAAGGAGGAAAAGGAAAGGCAGACGAATGGCATGGCTT TACTAAAGAAAAGATGTTGGCCTCATACTCTATACTCAGGGCTTAATGAACTGGAATCTGCATAACTC AGCAGTCAACCCAGAAGGGAATGGTTAACTGAGCTTGTATTGCCTCGGAGAGCCTAAGAGCACCCG CACACTTAATTCTACTCCTGTCTAGAAAAGCTGTGAGGAGTCGTTTGGAAATGCAATGTAGTTATTA AGGGCTGTTAACCAGCCTGCATTACATCTGGAAGTCAGGACTTGGGTGCTGACTATGAAGGGCCCTGTT TTCAAAATCTAACATTGCAAGTGTAATGGCAAGAAGCCTCCGTTGTGCTTTTTTTTCTCTTCAGT AACTTTTGCAACATTATTGCATAGAAGATCCCTGACCATTTACTAGGAACCTGGTTAAGCAAGCACTAA TCTCTTTTCTGGAGATCAAGGATGCAACCTCAGGTTGAGAAAGAAACAGGGTTCCCTGGGCCCATAG ACTGTTTGCAGGCATCACTGCTTCCCCCTGACACCTCACAACTAGCAAAAATTGCTTTGTCTTTGGA AATTATAGAGGGATTTGGGTATCCAGATTGTGCAGATGCAAACTTAGGCTGTCTTGATGCAAACTTAGA ACCACAGAAATGCTTTTAAATGCCTGTTTAAAGATGGAATTGTTGTTTTATAATTTGATTTTAGTGC TAAATAAATGATTGGCTTTGTACATGAATATGTTCTGTACAAGTCTTTCACTAGTACTACAGATAA TCAAAGCTATCAGAATTGTCTTTGATCATATTTGACGGTAATACACAAATCCATGTTTATGACGCT ACGCGTAAGCGGCCGCGGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001126130.2
Summary:	This gene encodes a protein subunit of the GINS heterotetrameric complex, which is essential for the initiation of DNA replication and replisome progression in eukaryotes. Alternatively spliced transcript variants encoding distinct isoforms have been described. [provided by RefSeq, Jul 2008]
Locus ID:	64785
MW:	56