

Product datasheet for **SC214609**

POLR1B (NM_019014) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	POLR1B (NM_019014) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	POLR1B
Synonyms:	A135; RPA2; RPA135; Rpo1-2; TCS4
ACCN:	NM_019014
Insert Size:	2000 bp



[View online »](#)

Insert Sequence: >SC214609 3'UTR clone of NM_019014
 The sequence shown below is from the reference sequence of NM_019014. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AACATCAAAGTGAACTGGATGTTGTTAACTTGATGTTGACCTTTTGGATTAAGAGGACTATCAGATT
AAAGCAAAATGTAATTTTAAATCAATGAAGATATCATTACCAGGTTACTCTTGAGATTTTCAACGGTG
TTAGAACTCTCAACCAAGACCTGAAAACCAAGTATGCAAGGTTTCTGAATCTCTCTGGTAGATTAACCTA
TTGACAATGATTTTCTGTATCTTTGTTCAAAAAGTTCATGTCTCTCAAAATAGAAATATTGATAAA
TGGAAGAGCATACGGTGACAAGTCTCCTTCCAACCCAGGTTCCCTACACCCTGCTCTCAGCAGGCAG
TGAGTGTACACACCTGTTAATCCATCTTGAGCAGGACAGTACTATACAAATAGAATGCAAGCTGTAAT
GTAATTTTATATTTTCTATAGCCAGTGAAGTAAAAACAAACAGGTACAGTGTTTTTTACCAGCTTT
ATAGAAGTACAGTTGTTACATATTTAATGAATACAATTTGATGGGTCTGACTATATGCACACACCTTTG
ATACCATCACCACAATCAGGGTAATAAACATACCTGTCATCTCCACAAGTTTCTCTGCCCCCTTTGTT
TTTTGCTTTTTGGTTGCTGTTGAGTTTTTGTCTTCTGTGGTAAGAACACTTAACTCAAGACCTA
CCCTCTTAAACAAATCTTAAAGTGACGATATAGTATTGTTAATCCAGGCACCATGTTGTACAACAGAT
CTTTAGACCTTACTTGTCTTGCATAACTGAAGCTTTATACCTGTTGAACAACCTCTCCATTTCCCTGGCC
CCTAGCAACCACCTTCTACCCTGTTTCTATGAGTTTGACTATTACAGATATCTCATATAGTGGGATCA
TGCAATATTTGCTGTGACTGGCTTATTTCACTTAGCATAGTGAATAAGATTCATCCATTTTGGGAG
CCAGGCATGGTGTGTCATCTATAGTCCCTGCTATTTGAGAGGCTGAGGTGGGAGGATCATTGAGTG
CAGGAGTTCAAGGACAGCCTGGTAATATAGGAAGACCCTGTCTTGAAGACCCTGACCTCAAGTGATCC
ACCCACCTCGGCCCTCCGAAAGTGCTAGGATTACAGGTGTGAGCCACTGTGCCTGGCCTCCGGTGAGTAT
TTATATTTAGTCTACACTTCCATACTTGGCTTTTTTCTGCTTTTATATTGATCTGCTTTCATAGCAGT
GTGTAGAGTGCCACTTATGTTTTCTTCTGTGTACAGTATTTTATTGTATGGATTTACCATCCCCTGT
GTATTTAAGTTGTTCCATTCTTTGGCCATTATAACTTTTTTCTGCAAAATATTCTGGTGACTTATCTTTG
GCCATTATAAAGTGTGATAATAGATCATCTTGTATATACTTCTGCAATTATAAGATGTTTTTGTGATGA
TGAAAGCTTTCTAAGGTAATTTCTTCTGAGTTTGGTTCAGATGGTCTGCTACTAATTATAAGGCAGA
TAGGAACAGACAGAAAAATCTATCTTTCAAGATAGGCTTAGCAGGATGGGCGCAGTGGCTCACGCCAC
TAATCCCAACATTTTAGGAGGCCTAGGCAGGAGCAATCACTTGTGCTGGGAGTTCTAGACTAGCCTGG
GGGAGACTTCATCTCTACAAAAAAGCAACAACGACAAAAAAATTAGCCAAGCATAGTGCCACACCCC
TGTAGTCCAGCTACTTGGGAGGCTGAGGTGGGAGGATTGCTTGAACCCAAGAGGTGCGAGGCTGCAGTG
AGCCAAGATTGTGCCACTGCACTCCAGCCTGGGTGACAGAGCAGGACCCTGTCTCTATTTTATAAATTA
AAAAAGGCTGGGTGTGGTGGCTCACACCATAATCCCAACACTTTGGCTCAGCAGATTGCTTGAACCCA
GGAATTCAGTCCAATCTGGGCAACATGGGAAACCCAGCTCTACAAAAAAATTAGCCTGGTGTGG
ACGCGT AAGCGGCCGCGGCATCTAGATTGGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_019014.6](#)

Summary: Eukaryotic RNA polymerase I (pol I) is responsible for the transcription of ribosomal RNA (rRNA) genes and production of rRNA, the primary component of ribosomes. Pol I is a multisubunit enzyme composed of 6 to 14 polypeptides, depending on the species. Most of the mass of the pol I complex derives from the 2 largest subunits, Rpa1 and Rpa2 in yeast. POLR1B is homologous to Rpa2 (Seither and Grummt, 1996 [PubMed 8921381]).[supplied by OMIM, Mar 2008]

Locus ID: 84172

MW: 74.7