

Product datasheet for **SC214639**

PABPC1L2A (NM_001012977) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PABPC1L2A (NM_001012977) Human 3' UTR Clone
Symbol:	PABPC1L2A
Synonyms:	RBM32A
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001012977
Insert Size:	1538 bp



[View online »](#)

Insert Sequence: >SC214639 3'UTR clone of NM_001012977
The sequence shown below is from the reference sequence of NM_001012977. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

[illegible]

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_001012977.3

Locus ID: 340529

MW: 58.3