

Product datasheet for **SC214647**

RGPD1 (NM_001024457) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	RGPD1
Synonyms:	RanBP2L2; ranBP2L6; RGP1; RGPD2
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001024457
Insert Size:	1459 bp



Insert Sequence:	>SC214647 3'UTR clone of NM_001024457 The sequence shown below is from the reference sequence of NM_001024457. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC CAAGAAACCTCCATACCCAAGAAAAAATAGAAAGCAACCATGTTCTACTATGGTGCCAAAAATAGTTGA ACATCGACATCAAAAAGATTCTTTCTGTCTGCTTATTACATATTCTGGTTGAATTAAGATGATGAA ATTATCCTAAATACAGTTGGCCCTCTGTATCCATGGGTTCCACATCCACAAATCAACCAACTGTGGAT GAAAAAGACTTGAAAAAGATGAGAGAAAAATAACAATATGACAATAAACATAATACAAATAAAAAACAAT ACAGTATAACAACGATTACATAGCATTACATTATATTAGGGATTGCACAAGAGAAGGTATGTAGGTTA TATGCAAACTACTGCACCATTTTACATAAGGGAAGTTGAGAAGCATGGATTTTGGTAGCCACAGGGGTCCT GGAACCAATCCCTCACAGACACAGACGACACTTTACAGTAGATGAACACAAAGATGAAAGGAAAAAGTC TGACCTAGGTTTGCAGGAGAGAGTGAAGTCCATTTTGCAGGTGATGCCATTTTGTGTTTGACAT CGTCCCTCTGTAGTTCTTTCCATTCCAGTCTTGCACTCTGAAAGATACACTGAAGGAAAGTCCACACA GTGGTCAAAGTCTTTCACAAGACACACGTGAAGGTCTGCACAGCACAGTCACATTGAGAAAAAGATCT CATGCACAGACCCCTGTTTCTGCTTTCTAAAAGATCATCTTTTGACCTGCAAAAAGGCTGCAGTAA ACTGGGCCATTCCATACATTTGATTCATGTATTCAAATGCTACTTATGAGCTCTCTGTGATTGAGCCCC ACCAGGATGCTGGGGACACACGGGGCACAGTGCAGCCAGAGGCTTGGCCATCATAGAGCTCACATTCTA GTGGAATCAGACAGGGGGTTACAGGAAATATTCAAGTAAACAATATAAAATGAGATCACTTGATGATA CATGGATACATGTTACATACATGATGAGACAGAGTAAGAGAATGGGAAAGGGTTTTAGATGGGGGAA TGGCAGGCTTTCCGAGGCTGAAAGAGATCCAACCCACTTCTATTTCTTCAGGCTCCAGTGATATGAA AAAGTGCCAAACAGGACTATAAAATTTGGTCTACTGTAAATGGTTTGCAGTAAAGTTCTCAACCAGC TTCCAATCAACCAACTCACCAGTTAGACAATGGGCCCTCTCCTGTGTTTAAACAGTCCGGTGCCCAACC TGGGTCACATGCAGCTGGCAGGTGTTGCTAGTGTGCCCAACCATGTCTGTTCCCCCAGTGAGTTA TGTGCTTACAACGAGTCGTTTCTTGCAATTTTCAACTTTCTTACATTATATTTAATTAAGATGATC TAAACTGAGAAATGAATGTGAAAAGAAACAAAGCTCTTTTTGTTTGGTTTTAACTAAATTAATGTT TTTAAACAAGT ACGCGTAAGCGGCCGCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001024457.4</u>
Locus ID:	400966
MW:	55.5