

Product datasheet for **SC214655**

ZNF764 (NM_033410) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF764
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_033410
Insert Size:	1460 bp



Insert Sequence:

>SC214655 3'UTR clone of NM_033410

The sequence shown below is from the reference sequence of NM_033410. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
TACCCGGAGATATTCCAGGAGTGTGGGTGACGGCCTTAAAAGTGACTATCTAGACATTACGGGAGGCC
GAGATGGGCTCAGGGGCCGAACCTCTGCGGCGCCTGCAGGGAGGTCCCAGAATCTAGGGCAAGAGCT
GGCCCTGGGGTGGCAGAGTCTGATCTTGGGCCCTTAGCAGCTCCTCTGCCAGCACCTTGCTCCCCAC
TGCCAGGGCTCTCCAAGGCCCTCTACTGAGGCAGGGCTGAGTTGAGAACCCAGACCTACATAAG
GGGAACCAAAAGCTGTTTCTCAGCCAGAGGTGCTAAGAGGACTGAGGTAGAGGAGAACCCTGTTTAC
CCATTGCTTTTTAAATTTTTATCATTATTTTATTCTTGACCTCCAGAACACAATACTCATTG
TCTTTTATTTATTTATTTATTTCTGAGACAGCTTCTGCTCTGTCAACCAGGCTGGAGTGAGT
GGTGTGAACAGGGCTAACTGCAGCCTCAACCTGCAGTGCTCAAGGAATCCTCTCCATCAGCCTTCCAA
ATATCTGAGACTACAGCGGGCAACACTACGGCCAGCTAATTTTTAAATTTTACATAGAGATGAGGTTT
TCACTGTTGTCCAGGCTGGTCTCAAACCTCTTGCTGAAGTGATCTGTGCTCTCGGCTTCCGCAAGTG
CTGGGATTACAGCGGTGAGCTACCACCATGCCTGGCTTCCATTGTCTTTGTGTGTGTGTGTGTGTG
GGGTTTGTCTTGTGCCCAGGCTGGAGTGCAATAGCGCAATCTCAGTCACTGCAACCTCTGCCTCC
AGAGTTCAAGCAATTCTCTTCTCAGCCTCCCAAGTAGCTGGGACTACAGGCACCCACCACTCCATT
GTCTTTTACTGCCAAGTTAAATCCTGTGACTGTGACTGGATCTGGCCTGCGGTTGTATTTAT
TGATCTGAGTCCAGGCAGTGTTTTAAGTGATATGAGTTAGTAGTAAAATGAAAAGGATATTTACA
ACCAATACAGGTTTTCTGGTGTCTTTGAAAAATCCAAGATTTGGCAACTCGACCTAGGCTCCGACA
AAACCACAGTTGTTGGGAGCTGAGCCGGAGCCAGAACTCACAATTCACCGTGGCCCGGCCACCCATCCC
CACCCTCCCGTGGGCTCCTGACCAAGCACCACGTGCCATTATCATGGCTTGTCTTTAGTGTGTTTGG
TGTTTATCTTAAAGTACAGCGAAGAAGAAAGGGAAGTATTTCTTAGACCATGTCTCTGTCCAAATAGG
CCTTTCCATCTGGCATTTCCTGCTCATTATATGACTTGCTGAGTTCCTATGGTCATTGAGTTTGCAA
ACAGTGGTCTGTAGCACTCAGGTTTTCTGGAGTTGGAGAGAAAGGTGGCAAGAGAGAAGTAAAGTTG
GTTTTATTTTA
ACGCGTAAGCGGCCGCGCATCTAGATTCTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_033410.4](#)

Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID:

92595

MW:

54.4