

## Product datasheet for **SC214658**

### SLAP2 (SLA2) (NM\_175077) Human 3' UTR Clone

#### Product data:

Product Type:	3' UTR Clones
Product Name:	SLAP2 (SLA2) (NM_175077) Human 3' UTR Clone
Symbol:	SLAP2
Synonyms:	C20orf156; MARS; SLAP-2; SLAP2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_175077
Insert Size:	1681 bp



[View online »](#)

**Insert Sequence:** >SC214658 3'UTR clone of NM\_175077  
 The sequence shown below is from the reference sequence of NM\_175077. The complete sequence of this clone may contain minor differences, such as SNPs.  
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AGAGCTGGACAGCTCCCTCCTGTTTTCTGAAGCTGCCACAGGGGAGGAGTCTCTTCTCAGTGAGGGTCT
CCGGGAGTCCCTCAGCTTCTACATCAGCCTGAATGACGAGGCTGTCTCTTTGGATGATGCCTAGGCCCA
AAGGAGAGGCCAAAAGGGAAACCAAGGCTGCACACCTAGAACCCCAATTCAGCCTCCTGGGCACCCAG
AGGCAAGGCTGTGCACTCAGGGAGGGAGGGTGGGACACAGAGGTGCATCTAGGGTCCCACCTGTACCCT
TGCTCTTTCTCTTAGCCCTTAGAAGTCACTACTTCTTCCAGTGCCATGATCCCACCTGCGACCT
CTAGTGCAGTGCAGAGAAGGTGGGACCAGGGCCAGGGTCCAAAAGAGAATAAGCCTCCTGGGGGT
CTGACCTAGTTAGTTCTTGAGTTTGGGGTTCCAGTACCATCTGGATGCCCTGCCTGTTGAGCCCCATT
CTACATCCCACCATTAAACCAGGCCCCACCACAAGGTAGAAACAACCCCTAGAGTCAACGAGAAAGTC
ATTTTCAGAAAATCTACAAGTCTCGTTGAGACCACCACCATACCTCAGAAAGGTAGGACTGTGGCCTAGA
AGGGAAAGGAAAGCTGAGATGATGTCTTACCGTAGCAGCAGATCTTGGATGGTCCAGGCTCTATGTGAC
CTCCAGAGCAAAGAGAAAGACTTCGGACAGTCTAGGTCCTCAAATGTCCCCCATTGAGGACAACAGCCC
CAGCTCTTTTTCTTTTTTTGAGACGGAGTCTTGCCCTGTTGCCCATGCTGGAGTGCAATGGCACGATC
TCAGCTCACTGCAACCTCCATCTCCTGGATTCAAACAATTCTCCTGCCTCAGCCTCCAGAATAGCTGGG
ATTACAGGCGTACACCACCATGCCTGGCTAATTTTTTTGTATTTTGTAGTACATGGGGTTTACCACA
TTGGCCAGGCTGGTGTGCAACTCCTGACCTCAGGTGATCCACCACCTTGGCCTCCAAAAGTGTGGGA
TTACAGGTGTGAGCCACGGCACCAGCCTAGCTCTCAGATCTCTATTTTCATTTTGTGGCTTACCATTCC
CTAGCACACTGGCCTTGCCATCTTGTGGCCGAATAAAAAATAACACCTCTTAAGTCTAGCACACTGCAG
TGAGGCCAGGCACCTCAGTGCTGGGCAGGGCATCAGAAGGTGCTAAGCCCTCTCTCCACAATGCCAAG
ACGGAGACCACAGCCTACACCAATCCAGCCCTTGATTTCCCTGCTGCCTCCATAAACAGAAAGAGGTC
TGCTGGATCCGCTAAGGGATCAGGGAGAGGAAGAAAGAGGGATGGGGTGGGAGGCACCCCTCCAGTGC
TCCTACTGGTTCCCAAGCTACAGGTGGGGTGGGAAAGGCTTTATCAGGTATCATCAACAGGTTCTCAAT
TAAAGATTTGATTTATTCAAGTATGTGAAAACATTCTACAATGGAACCTTATTAGATGCTGCATGTA
CTGTGCTATGGACCACGCACATACAGCCATGCTGTTTCAGAAGACTTGAAATGCCATTGATAGTTTAAA
AACTCTACACCCGATGGAGAATCGAGGAAGACAATTTAATGTTTCATCTGAATCCAGAGGTGCATCAAA
TTAAATGACAGCTCCACTTGGCAAA
ACGCGTAAAGCGGCCGCGCATCTAGATTGAAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

**Restriction Sites:** SgfI-MluI

**OTI Disclaimer:** Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:** The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:** [NM\\_175077.3](#)

**Summary:** This gene encodes a member of the SLAP family of adapter proteins. The encoded protein may play an important receptor-proximal role in downregulating T and B cell-mediated responses and inhibits antigen receptor-induced calcium mobilization. This protein interacts with Cas-Br-M (murine) ecotropic retroviral transforming sequence c. Two transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

**Locus ID:** 84174

**MW:** 63.4