

Product datasheet for **SC214659**

ZNF569 (NM_152484) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF569
Synonyms:	ZAPI; Zfp74; ZNF
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_152484
Insert Size:	1470 bp



Insert Sequence:

>SC214659 3'UTR clone of NM_152484

The sequence shown below is from the reference sequence of NM_152484. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GTTAGACACCAGAGAATTCATACTCATTAGAAACCTATGAATATTGTGAATATGGCAAGGCCATCTGA
AGGAATTAACACCTCATTGCACATTACATGATCACTTCCAGAGTAGAAAATATGAATGTGGGATAGCC
TTCTGAAAAAGCCACAATTTATGAAACATTAGAGAATTTCTCCAAGGTGACAAATTATATAATGAAAA
AGCTGTTACCAGAACTTTACGACAGACATCTTATTGATAATTTAATCAGCATTCTCATTGAAATCTG
AAACAAGGTATCTTCCATCAGCATAACAACACAACTTGAGCTCCTAGCCAATGTATTAAAGTAAGGAA
AAACAAGTATCAACGAAGAACTATAAAAAATACTTTTATACAAGAAATAATTAGAAAAATGTTGATAGA
AAACCTAAGAATAAACACAACAATGAGTTAGGTCTTTATCAAGTTATAAAATATGTTGAAAGGCACAAA
AGAACACCTGAAAAATGATATTTTGTCTGAAAAGACTTGACATTGAAAAATATAATTTCTCTCCAAA
TCACTAAGTTTAATATACTTCTGAAGTTTACATCAGTTTTATTTTAAAGGAACCTCACAACTTATCCAA
AAATTTTATGTGGAAGTAAGAAACACAAATAGCTGAAACAACTTGACAAAGCAAGTAGAGAAGAC
ACACTCTATTATTAGATATATCATGTACTTATTAAAGCTGCTGTCATTAAACAATGTATTACCAAGT
GAATAAATCAAGCAGTAGTAAGATAGTGAGCATTACATAAAAAATCCAAGAAGATACAGGAATTCGTG
TATGACAGTGCTATTGAAATCACTGGGGAATGAATGGCTTATTTAATAAATGGAAGGTTTCATAGGAGAA
AAATTTTAGCTTTCTACTTTATATCATATACAAAAATAAAATCCAATAAGACCAGAAATATAAAAAAGC
AAAAATTTAACTATGATGAAAGTAGAATCCGTTACTTTAGGGTAACAAGATTTCTTTAAAAAACATG
TAAACATGGAATAGTCTGTTGAATTTTGGTGATTGGCCTCACAGGCACCCACAAATGTACACAAGAG
ATATATATCACTGTGTTCTGAATCATGGTTTGTATTGCCCAAAATGGATATAACAAATTTTATCAA
TTGTGGAATAGATGAATAAGTGTGATCTTTCTGTGATAAAATACAGAATTTAGAAACAAATCTACCT
ATAGTTATATGGATAGATTACAATATATAATATTGAGAGAAAAAATTAATATGCACAAAGGAACACAGT
ATAACATTTGTATAATTTAATAAGATCTTGCAAGAAGCAGTACTGTGATTGGTCATGGATCCATATAT
GACTTGGTTGCCTTTTGGGATGTGGGAAAGGTTTGGTGGGAGTTGTATATTTATCTGTAATATTTCTT
TAATAAAAGAGTGAAAGCAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_152484.3](#)

Summary:

May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID:

148266

MW:

58.4