

## Product datasheet for **SC214687**

### **GABPB2 (GABPB1) (NM\_005254) Human 3' UTR Clone**

#### **Product data:**

Product Type:	3' UTR Clones
Product Name:	GABPB2 (GABPB1) (NM_005254) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	GABPB1
Synonyms:	BABPB2; E4TF1; E4TF1-47; E4TF1-53; E4TF1B; GABPB; GABPB-1; GABPB2; NRF2B1; NRF2B2
ACCN:	NM_005254
Insert Size:	2000 bp



[View online »](#)

**Insert Sequence:**

>SC214687 3'UTR clone of NM\_005254

The sequence shown below is from the reference sequence of NM\_005254. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CGTCTTCAGACTAATAAAGAAGCTGTTAATTGAAATGAACATGTAGTTTACTTTTGGTCAAG
AAAGAATACAATCTTGAAGTGTACACAACAAAGGTACAGCCATGGGAATACAGAATGATAGAAGAGACT
ACAGATGGATAATTGGACTTAAGCCATGAGCTCTGAGTCTTGTAAACATAAAAACCTTACTTTAGAAAGTT
GTGAAATGATTTAAAACGAATCTGTAATAGTTTTTTTTTTTTTACAGTTCCAAATGAGTTGATAA
AGATTGTTGAAGAGATCCAAAACAGAATAAGCCACTGTTTTTGTGAATCTTTTGTATTTAGTACAA
ACCTTAATTTCTCAGAAACGGAACAGTTTTAAGGGTGATCGTTGTTGGTTAGGCCAAATGTTGTGAAT
AATTATGGTGGACTGATGCTGGAATTACTCTGTAGGTATAAACCTCTGTATGAAGAGAAGATTTCTCC
CAGGAAATCTTTGTACAGCTTTAAGTTGTGTCAGATTCTCTGAAAACATTTTTTAGAAAGCAAAATTTT
TATATTTGTTCAATTTACGCTATACCCAAGTAGATTTACATGTATATGAAGCAAATTTTTTAAAAAATT
TCTGTTTGTACATATTCTGCATGTTTTATAATTTCAAAATGCATCACTTACATAGGTATTTCTCCACA
GAAATGATGAAAGTGACCAGAAAAAACAAAAACAAACCCTTACTCTGTAGGTCATTGAAACGAAG
TAAGCTGGCAGCTGGTTTTATTGGAATGACAGTGTCTCGGAAGGAGCAGCCTACAAGATAACTTGAAAT
TTGCCAATTCTGCAAAATCTGTGCTTTTTTGAATAATTAAGAGTGGGGACGTGAAACTGTATTCTGTGC
CTTCCATCATGATTTCCACATGAAAGCACTTTAAGGCACTGATTTAAGATAATGTTTTTGGAAAACCC
AATGCATATGGGTTTCTGAAATATTTTATGGACTTATTTCTCCCAGGAAATGATTCTTACGAAAAAAA
ATTGCTTTTGTATGTAGAACAGGAACTTTTTGTATTACAGTGTGCAATAGACATGCTAATGTAACCT
CTACTTTTCTTTTGAAGCTCAGTGTCTGTGCTATGACTTGCTCTCATCACAAATTTGTTGAATTTCCA
CAATGTATGGACATTAACACTGGCAGACTGTTCACTTTTTTTTTTTTTTTTTTGGTAAATATTACTT
CAAACCCCTTTTTCTTGCTTTATTTTTTCAAGTGTATTGCTTTTATGAACTGTTTAAACCCTGAAATCCC
TCTAGGTTATCTATACTGTATAAAAAAGCAATTACCCCTAAAACCTGACTCTGGCCTACTTTTCTATTT
TGCAATTAATATCTTTTTTACATATGTTTATTGACTTATGTTTTTATCACATCTTATTAACACAT
TAAAAATGTTATCTACTGCATATTTGAGCCACCATATGAAAATATTTTACTCTCCACAACATGACT
TAAAAATCTGAGTAACTCTTGAGAATTTGTTAACTCTTCACTGGAATTACCAGTATATTTGTTAA
TGTAATAAAGGGTCTATGGTCTGTCAGCCTGATGGTACTATTGGAAGCACTTAAAAAAGAGTTTAGG
GCAAAATGAAGAACTCTGTAAGGAAGAGGAAAAAGATCTTTTTCAATTTTCCACTAAATTTGACATATT
TTAAGTGTACAAAAACTAGTTTTTAAAAATGCCAGGCCATATTTTCAATAGCATGAAGTGAATTTTTC
TGCAGTCTGAAAAGTCAATAGAGAATGTTGTTTCTGGTACTTCGGTACATTAATTTGTGATGTTAGG
CAGGGTGTGGTACCACAGCCAGGTACCCTGACATGTGGTGTGGAGCAGAGAAGGGCCAGTGTCTTCA
CCCCAGAGGGCAATCACTGTAGAGTTTTGCTTTTCCAGACCAGGCTGGCCCTGAAAAGACTCAAGACAGT
ACGCGT AAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCACCGCCGCTTCTATGAAAGG
    
```

**Restriction Sites:**

Sgfl-Mlul

**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:**

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**RefSeq:**

[NM\\_005254.6](#)

**Summary:** This gene encodes the GA-binding protein transcription factor, beta subunit. This protein forms a tetrameric complex with the alpha subunit, and stimulates transcription of target genes. The encoded protein may be involved in activation of cytochrome oxidase expression and nuclear control of mitochondrial function. The crystal structure of a similar protein in mouse has been resolved as a ternary protein complex. Multiple transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

**Locus ID:** 2553

**MW:** 78.1