

Product datasheet for **SC214695**

CLIP170 (CLIP1) (NM_002956) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CLIP170 (CLIP1) (NM_002956) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	CLIP1
Synonyms:	CLIP; CLIP-170; CLIP170; CYLN1; RSN
ACCN:	NM_002956
Insert Size:	1436 bp



[View online »](#)

Insert Sequence: >SC214695 3'UTR clone of NM_002956
 The sequence shown below is from the reference sequence of NM_002956. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACCAACTGCAATGACGACGAAACCTTCATGAAGCCTCCAGTGGAGAAGTGGGCTTGCTCAGACGCAC
TCGCATTGACACAACTAACACCCAGCATTGTGTGTGCAGACTTCAGGAGAAGTCTATGTTATTTTTAAC
CCCGTCAACAAATCTAGGAAAATATTTTGTACTTCAACAAATTGCCCTTATGCTCCCGTATGAGTTA
GAATAATAAATATTTAGTAGGTGAGTTTACCTCGAATTTTGTCTTCTGATTTTACGTTTGAAGAC
ATTGCACCAGATGCCATTACATTTATTGGCCCCCGACCTTGTAGAAAAACCCCTACCCTCACAATACC
TTATTTAAGTAACTTTAAATTATGCCGTTACTTTTCATATTTGCACTAAGATATTTCCAGGCTGCATTT
GTATATTTAGATTTTTGGTTAAGCTTTGACACTGGAATGAGTTGAAAAATGTGCCATTTTGCATTTT
CATCTACTCATTTAAAGTATTTTATTCTTATTCAAAGAAATATCTGAGCTCTTTCGACTACCTGTTATC
AGTAGTGCTTTACTTCAGGCTTGATAACTTAGGTGTGATTATAAAATCATGAAGCAGGTAAAGGGA
GGGGCAAGCCCCAAACTGCTGTGGGGACATTTATAATCTATATGCTGCACCCACTTAATCTACTGTG
GTGTTTTGTTTATTAGTTTTGCATAATTTAGCTTCTATATTTGTATGTATATTTTTTAAAAATCT
ATATTTGGGAAAAAACATACACAATGTGTCTTTCTTTTGGACATTTACCTTTTTGAAAAAGAAAAAC
ACTTAAATGATCATTAGGACATAACAGACTAGGCCAGACATAGCATCTTGTGGCTTTGCAACCATTTT
CATTTGTTGTTTTCTTTTATTTCTTCCAGATTTAAATAACTGGAGGAATTTTCTCCAATTTTTTTT
TTCCTTCTGCGCAGGTATCCCCAGCAGTCAATTAACAATAAGCCAGTATAAAACACCTAAATAACCAA
TCTACAATCTCCCTTACAAGTTTTTTTACTGTTTTTATGATGAATGTACGATGAGAAATCAACGTTAA
TAATTCTGGATTTTCTTATCACAAAAAAGAAAATGAAGGACCTCAAAGCACCTGAACAGTTTATCGACC
AGTTTGAATCTATTTATCTTCAATTTGAATGTCTTCTAGAATATGTAAAAAGTCATAAAATGTATCTTCC
ATGCTACATGTACAATAAGAAGTCTATAATTGTATATATGCCTTTGATGTATTTTCCCCTCAAGATTA
TCAACTGTGTGTTGACAGTGAATATCAATCTGGTACCAGTTGAAATTTTTGGTTATAAATGTAATAC
GAATTGTTTCAAAACAGAAAACATGTAAGCAGTATTAATTTGGCCAAACAAG
ACGCGTAAAGCGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCACCCGCCCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_002956.3](#)

Summary: The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011]

Locus ID: 6249

MW: 55.7