

Product datasheet for **SC214699**

SNX24 (NM_014035) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SNX24 (NM_014035) Human 3' UTR Clone
Symbol:	SNX24
Synonyms:	PRO1284; SBB131
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014035
Insert Size:	1486 bp



[View online »](#)

Insert Sequence: >SC214699 3'UTR clone of NM_014035
The sequence shown below is from the reference sequence of NM_014035. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ATATTTTACCCTCATCTACAGCCAGGTAGAAATCCTACATGGCTAAAAGAAGCAGAAGCAAGTTTCGA
AGTCACAGTCAAGGAAATCAATACCTACCAATTTAACCTAAACGCTATGATATATAACAGCTCTAGCTA
GTGGTAAAGTGCACAGTCCCAGCTTAATTCAGGGCAGGGACATTTCATTAGAATGGTGCTCTTAAAAA
TAGAACTGAACCGGGCGGTGGTCAGGCTAAGGCCAAGTGTAAAGAAGTAGAGTGTAGCTGCCAGCG
TAGAAACCATGAAAAGGAGGCCACAGGAGATTCTGGGAGCACTGGGTGTAGCAAAACAAAGCCACTC
TCTGCTTCAGTCGCACCATTTGCTAATTGAAAATCATATCCTGAATCATACTGAGACTGATCAACTTTG
GTAGCTTTTTTTGTTGAGATCTTATGACACACTACTCTTCTCACCCTGAGATTTTCTCAGCCAGTGATAG
TACATTCTGAAATGCTGGCACCAGGAGACGGCCACAGACACACACTGCTAAATGTGAAGATGGAAGTAA
ACTGGAATTAATTTATACTGACAATATTATGGCATTTTTAAGATCATGGCATTTTAATTTACATTAGA
GTGGAGTTGCATCATACTCAGGGTTAGCTTCCAGGTCACTACATAGGTAAATGGGCTATTAGGATG
ATCCTTGAAAGCCCTTTAGAAGGGTGCCATGTTGGAAACCTGTACATCCACAACAAGTAGCTTTTCTC
CTATGTTGGAAAAAAGACTGTTTCTTTGTTTGAAGACCAAGTGAAGTGTGGTGTGTTGTTAGGGGC
CATTTTGTAAAAAAGACTGTTTCTTTGTTTGAAGACCAAGTGAAGTGTGGTGTGTTGTTAGGGGC
CCCCGCTAGTGGCAGATGTGAAGTGAAGGAGTGAAGCGCCACCCAGCGGATGGACAGCACCCACCT
GGGTTTACTCAAGGTGTGCATTCAATTTAGGTGGGATCGCCACAGGATTCATGTTATTTTCTTACG
GCTTCCTTTTCACTGACCTCATTTGTTGAGTTAATGTAAGTTAAATGTGTTTATGATATAACTCCACTG
TACATCATCCTTTGAGTAGTAAAGGATAAAAGCATATATACTACCTATATGTATGTGCTGTATGTGGGC
ATTTCAATGAGATCTAATTAATAGCTAGCTATTTATGTTTATTCGTTTATGTAAGTTCTGTGGGAGCA
AGGTATTTAAATCAAACTAATAACTACATCATGTTTGTGATTAGGATCTAAATATTCAGGTTTTAA
GCCTGCTGCAAACTTTTAAATATTATGATAGATTCTGTACTACATGTGGGAAACAAGCAAGAACTAAA
TAATCAAATGTTGTCAACCAAAAGTAATAGTTGGGTATTGGAGATTTTTTAAATGTTTTATGTTAT
TAGCTATTTGGAGTTAAATAAAACAGAACAGGAAA
ACGCGTAAGCGCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_014035.4](#)

Summary: May be involved in several stages of intracellular trafficking.[UniProtKB/Swiss-Prot Function]

Locus ID: 28966

MW: 57.3