

Product datasheet for **SC214709**

SCOC (NM_001153585) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SCOC (NM_001153585) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	SCOC
Synonyms:	HRIHFB2072; SCOCO; UNC-69
ACCN:	NM_001153585
Insert Size:	1492 bp



[View online »](#)

Insert Sequence: >SC214709 3'UTR clone of NM_001153585
 The sequence shown below is from the reference sequence of NM_001153585. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
ACAACTGACACAAAAAGCAAAAGAAAGTAAGGGATTGACACCCTTCTGTTTTATGGAATTGCTGCTGAT
CATTTTTTCTTTAAAAGTTGGATAGATTCCAAAAGTTACAGTACCTTTGTGGCTTCATTGAATTTAT
GAAGATAATGTCAGATGTAGACAAAAAACAACAATAACAGGAGACTTCCATAAGTTTGTGATTATGT
TAGTCTATGAAAACGTGCAAAATGTATTGTAGAGACTTTATGATTAGAATTGCATATTTATGAACTT
AAAGATGAATGTTTTATTGAATTTGTAGGTTTAGCACTGTCTTTTATTATAGGATTAGTAAGATATA
AGAAAAAACCACCGTGTGTGAAAAAGTGACCAAAATCATGTAATAATGCACAGCTTTATGTACCCT
GTCCACCATCTTGTGCCTCTTCCATTTGCCTCTTCTCTATTTCCCTTCCGCTAAGGAAAAAAT
TGGTGTACATTTGTAAGTAATTTAATAGTTAATCATCTCTGAGAGTAACCTGTATTTAATTGTT
GAAACTTAACCAAAATAAGATACTGTCTCAGCTAGGGCTTGTCTTTGTGATTTAGTGTTAAGATAGG
AATGCTAGTGTCTCTTAATTAATTGGAATAGATGGAGGCTAAAAATGAAGTTTTTCTTTGAACTG
AATTAACTTGGGAATATTTGTTGTTAAAAAAGTTCTTTTCCCAAAATAACTCATTTTGTATTATCTGA
AAATATAAATTTCTGGTCATGTGTATGTTAAAAATAGAAAATTTTGGGAAAAATGGAAATAGGGTGG
AAAGTACTCGGTAACAGTAGTAACCAATATTTTCACTCCAGATTTGTGTTTTCTCTGGCACAGAGTA
GATCTTTTGGGAAATATATGAAAGTGGATTAAGTTTGACTACCCTTATGTTAGCCACATCTGGATGA
GAACAGTTACAAGAGTTTGGTCTCTAAGTTGATTTGTACCCAGTGGTCAACTTCTGCAAAATCCGT
AATGGTGTATTAGTATTAGAATAGTGAATAAAATGGGAAAGTTATACATGTATACTTATTATCTTGCTC
AGTATTTTATCTCACTTGTCTAGAAATTTCTGTAAACCCTGCTACTGGGTTTGAAGAGTTTTAGTCAT
CCTTTAACAAATTTTAAAAATTTAGCTTCTAGATTCCATTTGGTAAGGAAATCAATATTGGAAGTATTG
CTAAAATCTTATAATATGAAAAGAGATCCACTAATGTAGCTTAAGGTTATTAGATTTGGGCTTTAATC
ATGGAATAATCTTATGTATTGGTGAAGAGTTGATGAATGACTTTAGCTGTGTAATATAAATAGTCA
AACTGCAACATTTTGCATCCCTTTTGTGACCTAATTTACAGACATTTAAATTGTGTTGCAGTCTGCT
TTGCCGTTTAATAAAAAGCTATTTAGAGGTTGAAAAA
ACGCGTAAGCGCGCGGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_001153585.1](#)

Summary: This gene encodes a short coiled-coiled domain-containing protein that localizes to the Golgi apparatus. The encoded protein interacts with ADP-ribosylation factor-like proteins. Pseudogenes of this gene are found on chromosomes 1 and 14. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2009]

Locus ID: 60592

MW: 57.3