

Product datasheet for **SC214716**

ZNF160 (NM_033288) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF160 (NM_033288) Human 3' UTR Clone
Symbol:	ZNF160
Synonyms:	F11; HKr18; HZF5; KR18
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_033288
Insert Size:	1493 bp



[View online »](#)

Insert Sequence: >SC214716 3'UTR clone of NM_033288
The sequence shown below is from the reference sequence of NM_033288. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ATGCATACCGGAGAGAAACCTTACAAATAGGTGTGGTGAGGTCATTAGGTACAATCACTCCTTTTACA
TCAGTTAATTTTATTCTTGACAGAATCCTTACAAATGTAGTGACAGTGGCCAATCCCTCATGAGTTGAA
GCATTAATAGATATGAGAGGCCATAAGCAAGAGACATCATGTAAACATATGTGGCAGAGGGTCTATCCA
GGCCTCGCAGGTTACTAGGCATCAAGATTTATATCTTTGATGAAACGAAACAAATGTAATATGCATCCT
GAGGCCATTACCCAGTGACCGATGGTAAGTGAGGATTCCTAGGAGGAATAACAGTCTCTGGTTTCCCTG
TTTGCCTTTGATATTATACACTGTAGAATACTCACAAGTCCAAATATGCTAAAAATTATATATTTTAA
CTCACATACGAAAAGGTTGCAGGATATTTGTAGGCAGTCAGTTACCTTCACCTTATGAAATGTTTCACT
GAGTTATTTGAGGTTTTTTGGAAAGCCTACTATTGCGTTTCAATGTGAACCTTGAATCTTATTGTGCA
TCCTTACACACCTTCCATGGTGCTTTCTTGAAAGATCATTGGGATGGAAGGATCATTGATTGGGTGAA
GATCATTGATTAGGTGAAGGATTATTTCTATCCAATTTGTGAAGAAGGAGGACTTTGCTTTTAAATTA
AGTATCATCTGAATTAGCATTGGGAGTGGCGAAAAACAATGTAAACTATGATGTCACTCACCATTCT
GATAATGTTTCAGGGTGCCTTTCTCTACAGGAGAGTACTGTGGCTTAGAGGAAAGAAATGGTCTATCA
ACTGAACATGAAATGGAGCAGGCCAAGACCTTAGGACATTGGGATTTTGTGGGAGGAGAGTAATAGGT
AATTAGACACTGATTGTGTGGTAGAAATACTGCAGGGGAAAGGTCGCCCTCTTATGCATCAAAGAGCA
ATACCTGTTGTTTAGCAAAGAGTGATGAAAAATTGATCTTGTTTTGAAATTGAAGAGAGAGGCCAGGCG
CGGTGGCTCACACCTGTAATCCAGCACCTTTGGGAGGCTGAGGCAGGTGGATCACCTGAGGTCGGGAGT
TCGAGACCAGCCTGACCACATGGAGAAACCCCAATTGTACTAAAAATACAAATTAGCCGGCGCTGGT
GGCAGGTGCCTGTAATCCAGCTACTTGGGAGGCTGAGGCAGGAGAGTGCTTGAACCCAGGAGGCGGAG
GTTGTGGTGAGCCAAGATCATGCCACTGCACTCCAGCCTGGGCAACAAGAACAAACTCCATCTAAAGA
AGAAAAAGAAAACTGAAGAGAGAGTAGTTGTCAGAGTGAAGAGTTTAAATCAACAGAGTAGAGTAACA
TATTCTTTAGTAGAACTCACATTATCTGCTTGAAGCACATTTTGCTAATGTTTTATTGAGATGTATCA
TAGATATAATGTTTTAACTAATAAAATAGAATGGTTTTTCTCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_033288.4](#)

Summary: The protein encoded by this gene is a Kruppel-related zinc finger protein which is characterized by the presence of an N-terminal repressor domain, the Kruppel-associated box (KRAB). The KRAB domain is a potent repressor of transcription; thus this protein may function in transcription regulation. Multiple transcript variants have been found for this gene. [provided by RefSeq, Apr 2016]

Locus ID: 90338

MW: 57.1