

Product datasheet for **SC214718**

ZNF160 (NM_001102603) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF160 (NM_001102603) Human 3' UTR Clone
Symbol:	ZNF160
Synonyms:	F11; HKr18; HZF5; KR18
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001102603
Insert Size:	1493 bp



[View online »](#)

Insert Sequence:	>SC214718 3'UTR clone of NM_001102603 The sequence shown below is from the reference sequence of NM_001102603. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC ATGCATACCGGAGAGAAACCTTACAAATGAGTGTGGTGAGGTCATTAGGTACAATCACTCCTTTTACA TCAGTTAATTTTATTCTTGACAGAACTTACAAATGTAGTGACAGTGGCCAATCCCTCATGAGTTGAA GCATTAATAGATATGAGAGGCCATAAGCAAGAGACATCATGTAAACATATGTGGCAGAGGGTCTATCCA GGCCTCGCAGGTTACTAGGCATCAAGATTTATATCTTTGATGAAACGAAACAAATGTAATATGCATCCT GAGGCCATTACCCAGTGACCGATGGTAAGTGAGGATTCCTAGGAGGAATAACAGTCTCTGGTTTCCTG TTTGCCTTTGATATTATACACTGTAGAATACTCACAAGTCCAAATATGCTAAAAATTATATTTTTTAA CTCACATACGAAAAGGTTGCAGGATATTTGTAGGCAGTCAGTTACCTTCACCTTATGAAATGTTTCACT GAGTTATTTGAGGTTTTTTGGAAAGCCTACTATTGCGTTTCAATGTGAACCTTGAATCTTATTGTGCA TCCTTACACACCTTCCATGGTGCTTTCTTGAAAAGATCATTGGGATGGAAGGATCATTGATTGGGTGAA GATCATTGATTAGGTGAAGGATTATTTCTATCCAAATTTGTGAAGAAGGAGGACTTTGCTTTTAAATTA AGTATCATCTGAATTAGCATTGTTGGGAGTGGCGAAAAACAATGTAAACTATGATGTCACTCACCATTCT GATAATGTTTCAGGGTGCCTTTCTCTACCAGGAGAGTACTGTGGCTTAGAGGAAAGAAATGGTCTATCA ACTGAACATGAAATGGAGCAGGCCAAGACCTTAGGACATTGGGATTTTGTGGGAGGAGAGTAATAGGT AATTAGACACTGATTGTGTGGTAGAAATCTGCAGGGGAAAGGTCGCCCTCTTATGCATCAAAGAGCA ATACCTGTTGTTTAGCAAAGAGTGATGAAAAATTGATCTTGTGTTTGAATTTGAAGAGAGAGGCCAGGCG CGGTGGCTCACACCTGTAATCCAGCACCTTTGGGAGGCTGAGGCAGGTGGATCACCTGAGGTCGGGAGT TCGAGACCAGCCTGACCACATGGAGAAACCCCAATTGTACTAAAAATACAAATTAGCCGGCGCTGGT GGCAGGTGCCTGTAATCCAGCTACTTGGGAGGCTGAGGCAGGAGAGTGCTTGAACCCAGGAGGCGGAG GTTGTGGTGAGCCAAGATCATGCCACTGCACTCCAGCCTGGGCAACAAGAACTCCATCTAAAGA AGAAAAAGAAAACTGAAGAGAGAGTAGTTGTCAGAGTGAAGAGTTTAAATCAACAGAAAGTAGAGTAA TATTCTTTAGTAGAACTCACATTATCTGCTTGAAGCACATTTTGCTAATGTTTTATTGAGATGTATCA TAGATATAATGTTTTAACTAATAAAATAGAATGGTTTTTCTCTA ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001102603.2</u>
Summary:	The protein encoded by this gene is a Kruppel-related zinc finger protein which is characterized by the presence of an N-terminal repressor domain, the Kruppel-associated box (KRAB). The KRAB domain is a potent repressor of transcription; thus this protein may function in transcription regulation. Multiple transcript variants have been found for this gene. [provided by RefSeq, Apr 2016]
Locus ID:	90338

MW: 57.1