

Product datasheet for **SC214740**

PABPC1L2A (PABPC1L2B) (NM_001042506) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	PABPC1L2A
Synonyms:	RBM32B
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PS100062)
ACCN:	NM_001042506
Insert Size:	1891 bp



Insert Sequence:

>SC214740 3'UTR clone of NM_001042506

The sequence shown below is from the reference sequence of NM_001042506. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

[illegible]

Restriction Sites:

Sgfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_001042506.2

Locus ID: 645974

MW: 72.2