

Product datasheet for **SC214749**

Deoxyribonuclease I like 1 (DNASE1L1) (NM_001009932) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	Deoxyribonuclease I like 1 (DNASE1L1) (NM_001009932) Human 3' UTR Clone
Symbol:	Deoxyribonuclease I like 1
Synonyms:	DNAS1L1; DNASEX; DNL1L; G4.8; XIB
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001009932
Insert Size:	1485 bp



[View online »](#)

Insert Sequence:	>SC214749 3' UTR clone of NM_001009932 The sequence shown below is from the reference sequence of NM_001009932. The complete sequence of this clone may contain minor differences, such as SNPs. Red=Cloning site Blue=Stop Codon CAATTGGCAGAGCTCAGAATTCAACGGATCGC GTTGCTGCTATCACTCCTGTCCCTCAGCTGTGCCCTGCTGCCTGACGTCCCCCTACCCCCCAGGGCC TGCTGCCTTTTGGGACTTAAACCCAGCCTCCCCGTCCATCCAGCCCTGGGGCTGGGGGCTTCAACTA TAGTTGCCCTGTGACTGTAGTCCACCCCTGCCTGCCTTGTGATTGGCTCTTGTCTTTGGTTGGGCT TGTGCCTAGATTAGGAGAGGAAGCCAGGGGCCCTGCACTCATGCCACCTGCCAGGTAGTGTAGTATCAGG AGTGGAGACAAAGTGGGCTCTGGGTTGGGGTAGGGGAAGGGAGGGTTCAGAAAGAGGAATGAAGATGTTG TATGACAAGAAGGAAAGTTACTGAGAACAAAAACCCAGATTGGTGAGATAGGACACTTGTGCAGCAGATA TGCCAATGGGCCATGTTTATTGTGGATTGGTAAGAATCACCAGGAAACCATTAAAGCCCCAATAGTACAA GGAGGGTGGTTAATCTGCTATATCAAACCTCCTCCCTGAAACCAGCAAAACACCGGAAACATTTTGGCTC ATTATAATCCGGTGAACAATGCAGTCAGGCCTGTTATAACCGCTGAGCAGCCACACTCGCACCTCCTGGG TGCTGTAGTCTGTGTTGGTACAGGCTTCTGCATGCCTGGTAAAGTCCAGCCAAGGCTGGTCAAGGCAACA TCTCCACACAGAAAATCTGCACCAAGTTATGTAAGCTAAAAAGCTGTGTGAACCCAGGTGTCCCGAAAGG GGCTGCAGGACACAGCAAAATGCCAGCAGCGTGCCGACCCCTCCCTTCCATCCTCTCCTCCAAAGAACA GAGGTGAGGAAAAACTGGCTGGGACGCTAGAAGGGTCATGTGTTAACTATAATCACATTTATGGTTTG GAACCATCACCCCAAGGTAAGAAAAAATAAAGGTATTCCCAGGTATGTTTGGCAAAATAAAATAAAGG TAATTAACCACTAAGAAATGGCTGTAACTGAGATTACAAAAACAGAGGCCTTTTGGGAGCCACAGG CTGGGGGCTGGGGCCAAACATACCCAGCTCTGGGACCAAAAGCCTTTGGAGTCCTGTCCCTGGATACAGA TGCTCAGGCCACCCACCTAGCACTCCTCTCATTTTTGCTGCTCCCCCAAAAGCCAGCGCTGTTCTACAAT AGTTACAGGTGAGACATCTAAGATCCCCAGGCCAAACCATGCCACACAGTGTCTGTTCCGAGGTCTTGTT ACAAATGACCAAGGTTAGCCTAGACCTTTGGTCAGGTGGTGAAGGCCCTGGGAATGCAGGTTCTAGAAG GGGTGGGGCATGGGGCCAGGAGGTTGCTAGAGATAGGAAATGGGGCAAATCGCAGGTGACCTCAGCAG TACCTGGATGATTCCACCAAGGAGGAATGGCCAAACAGATGGGTGAGAACCAAGCCCAATGTGATGCAC ACCTACCCTTGCTAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001009932.1</u>
Summary:	This gene encodes a deoxyribonuclease protein that shows high sequence similarity to DNase I. The encoded protein is localized to the endoplasmic reticulum and modified by N-linked glycosylation. Alternate transcriptional splice variants encoding the same protein have been observed. [provided by RefSeq, Jan 2015]
Locus ID:	1774