

Product datasheet for **SC214751**

RBM19 (NM_001146699) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	RBM19 (NM_001146699) Human 3' UTR Clone
Symbol:	RBM19
Synonyms:	Mrd1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001146699
Insert Size:	1492 bp



[View online »](#)

Insert Sequence:	>SC214751 3'UTR clone of NM_001146699 The sequence shown below is from the reference sequence of NM_001146699. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GACAGCGAGGAGCAGACCCTTCAGCTGTAGCTGGCACCGAGAGGGCTGGAGTGAGTGCGGTGATCTG GGCTCACCTCCCTGGTTCAAGCGATTGTCCCGCTCAGCCTCCCAAGGCTGGGATTACAGGTATGCACC ACCACGCCCGCTAATTTTGTATTTTTCATAGAGACTGGGTTTACCATGTAGACCAGGCTGGTCTCGA ATTCTGACCTCAAGTGATCCACCTACCTCGGCCTCCCAAAGTGTTGGGATTACGGGTGTGAGCCACCA CGCCCGGCTAATTTTGTATTTTGTAGTAGACAGGGTTTACCATGTGGGCCAGGCTGGTCTTGAATC CCAACCTCAAGTGATCCACCCACCTTGTCTCCCAAAGTGCTGAGATTACAGGTGTGAGCCACCATGCC CAGCCTGTTCTAACCCTTTCAACTGGGGGTCTCAAGTGGGTGAGGACTCCATGGCCACGGCAGCAGA ACTGTCTCTTCTGAAAACCAAGTCCGGGGCCCTGGGTGAGCAGCTCTAGGTCTTCCACAGACTTAC ACAGTTTAAAGAAAGAGCCAGCGAACATGGGGTGATCCTGGGGTGCCACTGGGATCCCAAGCCAGGCCC GGAGGTCTGCCTGTTTCGTCCCCAGAACTTGAGCTGGCATCCTCCGCTGGTTTGCACTGGGCACGGGG ACTGGAGAGCCACCAGGCCACTGAGCGCAGCATGCTATCACTTCTCAGGAAGTTGTGAGGGTGAGGAGA ATAGCCGGGCCACTCTGCAGAGGGGAGGAACAGGCACAGGCTAGAGGCACGGGTGATAGAGTGAGTC CTCTCCCTCATTATTTCATTTCATTTCATCCATCAGACACCTGCTGGGCATCTGCTGTGTGCG AGACACTGTGCTGGTGCCCCAGACACAAGACAGAACTCGGGTGCCTCTGTGTGCCTCCCACTGGCC TTCCCTGCTGCGTGGACGCAGGCTCTGTCTCTAGGCTTTAAGCCCCATGCATCAGGTGGAAGACAATT TCCAGCAGCCTCTGCCGCTTGTCT GTGACGGGCGCTCTCTCGCAGCGCCGCTGTGTCTCCGACGCTAGGAAAGGATTAGTGTGAGCAGA GCACGCAGGCTTTGGTGAACAGAGCCGAGCTTGGTTTCAACATCTAATTCAATTTGTTCTGCCGAGATT GAAATATGAAATGGTTGCCAGTCAAACGTTCAAAGGACGGGGAAGGGAGAGTGAGCGCATGCTGAAATT TAAATTTTAAATGAAATGACTTTAACCTTTGGTACACCTTGCAAAATGTAACCCATCTCCCATCTG TGAAAGCGGCTCTCGGTTAGAGCTTTCCACCCACTTCTCTTCCACTGAGCTTCACTAATAATATGT AATCATTTTTTTCAGCCGTCTTAATAATAAAAGATTAATAGAA ACGCGTAAGCGCGCGGCATCTAGATTGCAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001146699.2
Summary:	This gene encodes a nucleolar protein that contains six RNA-binding motifs. The encoded protein may be involved in regulating ribosome biogenesis. Multiple alternatively spliced variants, encoding the same protein, have been identified.[provided by RefSeq, Apr 2009]
Locus ID:	9904
MW:	54.7