

Product datasheet for **SC214780**

TRANK1 (NM_014831) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	TRANK1
Synonyms:	LBA1
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_014831
Insert Size:	1486 bp



Insert Sequence:

>SC214780 3'UTR clone of NM_014831

The sequence shown below is from the reference sequence of NM_014831. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
CGGAAATGTGGAAAGCAGAGAAAATACAAATGTCCACACAGCTGCAGCCTCCTCATCCTTCGGAACATT
CCATTCTGACTTAGAATTCTGAGCGCTGGGGCAGAAGACAAAAGAATGTGAATTAGCATTTTAAAAAT
AGGGGAGTCTAACAAACACCATTTGCTTTTTATTCTCATTACTTCTCTTCAGTGGCTTAGGATCTGATTG
CTCCCTCACCTCTAAAGTGTTCAGTAGAGTCTGGCAGTATTAGTTTCCCCACAACAAATTCAGG
TAGAGAAGTTGTAGTCGAAGGGAATGTTGGGACTCCTCACTGCATGGTCCCTGGGGTGTCTGGCTCATC
CCCATCATGACATGAGTTGGTAGCAGAGGACTAGGTGTTAGGACAGACTCTCCTCACCTCAGTTGCTGC
TGCTCTGGGCCAGGAAATCTCTATGGATTCCCTCATCTCCTGCCACACTGCTGCAGCATCAAGTGTG
CTATGAGAGCTGTATGTGCCAGGGGTGGGCCCTGCCACTGGCCACCCGAGCCAGGACGATTCTGG
TCATCTCTGAGCCCACTTGTAGTGTAGCTTCTTCTGTTAGAGATGTTAGCCTTTTCTCTCTGCGGAC
TTTTCTGCCTTCCCTTCTCTGTTTTTTGTGATCATTCTTAGCTACCTGTGCCAAGTCACAGTTAGACT
ACCTCAGAATGTAGTCAAGAATTCACCCAGTGTCTTTGGATGTCCCTTCCCTAGGGGGAAGGCTGGG
GAATGTTGCTGGGAGCATTCCAGCTGCCTCTCCAAAGTGGTCTAACTCTGCCTTGGCTATGGTGCAG
ACATTCAAAAAGTGTGTGTATACAGCCCTGTACATCTAGTTTTCTGGTTCTCCCTGGCAAGTTGCG
TACCCTAGCTTAGGCTGACGTTCTGATGAGTCTCTGCTTGCTGAGCTCATCACTGTGGGATGGTACCT
GCAGAGCTGCAGATTGGAGCAAGTGGACAGAGGAGGCCAGCTAAGCCAAATCCGTCAGTCCCACTT
GGCTGGCAGGCTTATAAGAGGGCGGATCTGGGCCGTTTGACAGGAGGCCATTGCCTACCTCTCCCA
CAGTGAGCCAGCAGCCATCAGCTGATCCACCTGGGAAACCTTCATGCCTCTCTGATGGTTACTGCCCA
CCCTTACCCACCCCTCAGCTCAGCCTGGTATGGAAAGCAAGGTGCACGTTGGTCTTTGATTGTTCTGT
CCTCACAGCAGAGCCAGCTTTCCCCCATGTTGCTGCTGCTTTGCTGCTCATGCTGCTTGTGTTTTT
TCTTCTTATCTTTCCAAATGTTTCTTGAATGCTGTTTCCAAATTTTCTTAGATTATGTTTCTAA
TTTTTATTCAGGTGTTTTTTTTTCTTACTTGACAGAAAAAAGAAACCCCAACAATATCTTGCTGTA
TTATCACATTTTTTATAAGTTAAAGCATTCTCTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
    
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_014831.3](#)

Locus ID:

9881

MW:

54.7