

Product datasheet for **SC214781**

ZNF213 (NM_001134655) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF213
Synonyms:	CR53; ZKSCAN21; ZSCAN53
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001134655
Insert Size:	1487 bp



Insert Sequence:

>SC214781 3'UTR clone of NM_001134655

The sequence shown below is from the reference sequence of NM_001134655. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CACGCCAAGATGGCCAGCCCGTGGGGTGAAGCAGCTGGCTTGGCCGAAACCCGGGGGAGGCCAGCCA
CGGCACATCCTGCTTTGTTCACCACTGGGACTCTCTTCCATCTGTGGCCACCTCCCGGGGTGTCCGAG
GGACCCAGGGTACCTCACACTCGGAGCTCGCCTGCCTGCTTGGCTCTGAGGACCTGCCAGCGCTCA
AAGGGAACGGAAGCCTTCCCTCCCGCCCCGATCTTGTCTCTTCCCTTCTGCGCTTAGCGTTCC
TCTTCCCTCTAGTTTCTGGAGCCCAACACATTCTGGCAGGGACAGCAGGGTGGCAAGGACTCAGG
TCTAGGTCCCTTCCAGAAGCCCCGAGCCTCATTTGACTGTGTGGCTCTTTGGCCCCACCTGTGGG
GTGGGTCCATGGGTGAGCCTCTGCCCTACCAACCTGTGCCTTTCAGTGGGCGTGGAGGACTGGCCTTG
GCCCCCAGGGGGTGTGGACTTTGGGAGAGACAGCCACACCTGTGGGACCGCGGTCTTAGTCACG
GCGGCAGGGGCTTCTGCCCCCTCCCACTCCCGTTTCCAGGCCATGACCACTCTGCCCTGTCTGGCC
ATACGGAAGTCCGCTGCTTTGCCCTCGGCCTACTTGCCCTAGCATGAGGCTCTGAGAGCCACCTGCC
ACCAATCTGGTGAGGATAATGGTGGCTCCAGCGACAGGAGGCCAACCTGGAGACCAAGAAGAGGGCGC
CTGGCTGCCATCTTTCTCCAGAGGTGGGGTGCACCACTAGCACTAGCACTCCATCAGCACTAG
CACCTCACTCCATCAGCACTAGCACTCACTCCATCGGCCCGGCCACCTGCTCCATCGGCACTGGCGC
CCTGCTCCATCGGCACTAATGCTCCACTCGGCGCCCCACTCCATCGGCCCGGCTCCATCGGCACTAATG
CCCCACTCGGCGCCCCACTCCATCAGCACTAATGCTCCACTCCATTGGCACTAACGCCCCAAGTCCAGC
GGCACTAATGACCGCTCCTTTGACATTGGTGGCCCACTCCATCAGCACTAACGCCCCGCTCCATCGGC
ACTGGTGTCCCACTCCATTGTCACTAACGTCCGGCTCCATCGGCACTACCAACCGCTCCATCATCACT
ATGTCCAGCTCCGTCGGCACTACCAACCTGCTCCATCATCACTACGTCCAGCTCCAACGGCACTGGTGC
CCCATTCATCGGCACTAACGCCCCGCTCCACCGGCACCAAGTGCCTCGCTCCATTGGCACCAACGCCCA
GCTCCACCGGTACTGGCTCCCTGCTCCATCGGCACTAACGCCCTGCTTATTGGCACTTTGCTGCTGCC
TCCTGAGCACTGCCTTCCATGAACAGGGACAGACCAGAGGCCCTGCAAGGACTCCCCCTCAGACCTCCA
AAGGGCAACAGAAGAGTATTAATAAACGTGAAAACTTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
    
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

[NM_001134655.2](#)

Summary:

C2H2 zinc finger proteins, such as ZNF213, have bipartite structures in which one domain binds DNA or RNA and the other modulates target gene expression.[supplied by OMIM, Apr 2004]

Locus ID:

7760

MW:

52.4