

Product datasheet for **SC214782**

ZNF213 (NM_004220) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF213 (NM_004220) Human 3' UTR Clone
Symbol:	ZNF213
Synonyms:	CR53; ZKSCAN21; ZSCAN53
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_004220
Insert Size:	1487 bp



[View online »](#)

Insert Sequence:	>SC214782 3'UTR clone of NM_004220 The sequence shown below is from the reference sequence of NM_004220. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CACGCCAAGATGGCCAGCCCGTGGGGTGAAGCAGCTGGCTTGCCCGAAACCCGGGGAGGCCAGCCA CGGCACATCCTGCTTTGTTACCACTGGGACTCTCCTTCCATCTGTGGCCACCTCCCGGGGTGTCCGAG GGACCCAGGGTACCTCACACTCGGAGCTCGCCTGCCCTGCTTGCTCTGAGGACCTGCCAGCGCTCA AAGGGAACGGAAGCCTTCCCTCCCGCCCCGATCTTGTCTCTTCCCTTCTGCGCTAGCGTTCC TCTTCCCTCTAGTTTCTGGAGCCCCAACACATTCTGGCAGGACAGCAGGGTGGCAAGGACTCAGG TCTAGGTCCTTCCAGAAGCCCCGAGCCTCATTTGACTGTGTGGCTCTTGGCCCCACCTGTGGG GTGGGTCCATGGGTGAGGCTCTGCCCTACCACTGTGCCTTTCAGTGGGCGTGGAGGACTGGCCTTG GCCCCCAGGGGGCTGCTGGACTTTGGGAGAGACAGCCACACCTGTGGGACCGCGGGTCTTAGTCAG GCGGCAGGGGCTTCTGGCCCCCTCCACTCCCGTTTCCAGGCCATGACCACTCTGCCCTGTCTGGCC ATACGGACTCGGCCTGCCTTTGCCCTCGGCCTACTTGCCCTAGCATGAGGCTCTGAGAGCCACTGCCC ACCAATCTGGTGAGGATAATGGTGGCTCCAGCGACAGGAGGCCAACCTGGAGACCAAGAAGAGGGCGC CTGGGTGCCATCTTTTCTCCAGAGGTGGGGTGCACCACTAGCACTAGCACTCCATCAGCACTAG CACCTCACTCCATCAGCACTAGCACTCACTCCATCGGCCCGGCACCCTGCTCCATCGGCACTGGCGC CCTGCTCCATCGGCACTAATGCTCACTCGGCGCCCCACTCCATCGGCCCGCTCCATCGGCACTAATG CCCCACTCGGCGCCCCACTCCATCAGCACTAATGCTCACTCCATTGGCACTAACGCCCCAACTCCAGC GGCACTAATGACCCGCTCCTTTGACATTGGTGGCCCCACTCCATCAGCACTAACGCCCTGCTCCATCGGC ACTGGTGTCCCACTCCATTGTCACTAACGTCCGGCTCCATCGGCACTACCACCCCGCTCCATCACTACT ATGTCCAGCTCCGTGGCACTACCACCTGCTCCATCATCACTACGTCCAGCTCCACGGCACTGGTGC CCCATTCCATCGGCACTAACGCCCCGCTCCACCGGCACCACTGCTCGCTCCATTGGCACCAACGCCCA GCTCCACCGGTACTGGTCCCTGCTCCATCGGCACTAACGCCCTGCTTATTGGCACTTTGCTGTGCC TCCTGAGCACTGCCTTCCATGAACAGGGACAGACCAGAGGCCCTGCAAGGACTCCCCCTCAGACCTCCA AAGGGCAACAGAAGAGTATTAATAAACGTGAAAACTTA ACGCGTAAGCGGCCGCGCATCTAGATTGCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_004220.3
Summary:	C2H2 zinc finger proteins, such as ZNF213, have bipartite structures in which one domain binds DNA or RNA and the other modulates target gene expression.[supplied by OMIM, Apr 2004]
Locus ID:	7760
MW:	52.4