

Product datasheet for **SC214787**

C9orf40 (NM_017998) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	C9orf40 (NM_017998) Human 3' UTR Clone
Symbol:	C9orf40
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_017998
Insert Size:	1495 bp



[View online »](#)

Insert Sequence: >SC214787 3'UTR clone of NM_017998
The sequence shown below is from the reference sequence of NM_017998. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAACGATCGCC
GAAGGGGCTGAGGTTGACATGGAGTCCATGTAAGGAGCCGAAGCAGTGGGATTGGCTGATTTGAGGA
GATGTCTCTAAGTGAATTCTCGTATTCTTAAGGAAAAAGTTATTTTCCATACTTGAAGTTATATTTCCA
AACCTGAGAAATGAAGAAAGATTGTTCTGACATTAATACTACAGTTACTACTGAACCTCTTAATAAG
GATTTGTCAAGGATAGAGTACAGTTGTAGGGGAAGTATTTATGTATGCATTCTTAGAGCAAAAAGTTT
TGTTTAAATTCTAGAATTGAAGTACTGATCTTATAAAAAAGAAATCTAGCAGTTTGAAGATAGGTGG
GAAAACTCAATATTCCTCTATCTGCACCAAAAAGTTATTTGTGGTATATAAAATGAATATTGTTT
TATAATAACTTGTTAATAAAGTACTTTCTAATACATTCTATTGACTCTGTTAGTTGAACAAATAGCTGA
CTTGAACATCTATGCAAACTTAAGATGGCGGGATTGTTGTAAGCTATTGTTTTAAAGAGCTTTCT
AAATGTAAAGTAGTGATAATTTCAATTTGGGTAGCGTGTTTGCAAGCTTCCAATATTTGATGTTGGTT
AAGCTCTACTATGGGCAACTGAAGATGGATAAAGAAAAATGAAAACTGAATCGGTGCCTGCTTCCCCCT
GTTTTCCAGGATTAGAGGAAAAATTTATTGTATAATCAGCTTCTTGGTTTTGAATTGCTTCGAGGCA
TGGTTTTATTCTTATTACTTTAGACCTGTAGTTTTCAACACTGACAGCACTTTAAAAATCTTTGCCTG
GGCCTCACTCTTGAGAGATTCTCATTTAATAGTTCTGAGGTGGGTCTGGATATAACTATTTTTTAAA
CACCTGTCCTGTTTCCCTCCATTCCTCTCATGTGCAGACAGGGTTGAGAACCAGTAGACTAATGGTCGT
TTTTCTGTTTAAAGGAGATAACTAATTTGAGCTGAAGCAATGCTTCTTAATTAGCTTTGTTTTGTTT
TGCTCTGTTGGTGGCTTTGTTACAACCTGAATTATTGTGTTATTACTATTTTATTGTTAAAGAAATAAG
TAAGCAATTTGTGATGTGAGTATCAGTGATTAAGTTAACTAACTTTTGTACTGCATCCAGAATGTTGGT
TTTGCAATTGAGTAACTGTTCTTGCTTGCAATTTTTTGTGTTGATGACATTAGATCCAAAAATCAAGA
CAAATGGTAAATGCCATTGAGAGGGAAAGAGAAAACTTGATTTTTTTTGTGTAATGAAGGATTTAAGA
ATGGGTTGACATTAATAAGAATGCTTTAGAACAGAAAGACAACTGTATTGCATTGTGGTCAGACATGGT
TCAAAGTCTTGACTGCCACTTCTACCTATGTATCTTTAAGCCAGTTATTTTTCATCTCCAAGCCTCA
AATTTCTCACCTGTAATAAGAGAAATAATAAATAGTATCTACCTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_017998.3](#)

Locus ID: 55071

MW: 58