

Product datasheet for **SC214805**

ZNF146 (NM_001099638) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	ZNF146
Synonyms:	OZF
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_001099638
Insert Size:	1479 bp



Insert Sequence:

>SC214805 3'UTR clone of NM_001099638

The sequence shown below is from the reference sequence of NM_001099638. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
ATTAGACACCAGAAAATTCATACTCACAATAAACCCCATGAAAGCCTTGAAAGTGGGAAAGCTTTTATT
AGAAATTTGCACCTCATCATGCCCCAGAAAATAATCCTTCTGAAGCAAGCACCACGAATGAGGTTAACT
TTAACAAGTACTAAAACTTAAGGGACACCAGAAAATTTGTTACTGAAGAGAAAGACATGCATATGATTA
AAACCTGTGTCCAACAGAGAAACCTGCAGCAGAGATAATGGTGAAAGTTTAGGCACATTTTCACTAAA
AGTGGGAACAGAAAATGGAGGCTGTTTTTATTGCTACCACCATAGATCTGGAGATCTTCGCCAGTAA
CAAGAAAAATTAAGTTGTAATATTGGAAGGAACAGACAAAAATAGATGACATGGTCATCTACAATA
ATATTAAGACTCCTGTGGTATTGATTGAGCAGAGCAGAAGAGATTTCAAAAAAGACCTATGCATTTAT
TAGAATTTGGAATATGATACAAGTGGCATCAGGAAATGAGGAAATATGGAACATTTTTTTTTAAGTGG
AGGTAGTTTGTCTCCAAGGGGAAAAATAGACCATGAACATACACAAAAGTGAATCCAGAGATTAATA
AACAGAAATATGGAAGCTATGTACTAAAATGCATATGACCTTGGGAATAGGGAACATTCCTTGGAA
CAATTTGAGAAATCTTAGCGGTTCTACTATGAGTATGAATGCTCAGGAAATACGCACTAGGATATTT
ACTGTGGCTTGATTGTAGTAGCCAAAAATAGAAACAATGAGAAAGCCCATCACGAACAATATGGGA
AATTCCTTATGTAATACTATATCTGTAAATGCAACAAATGAAATCCACATGTATTAAACAGAAAGATC
CCAGAATTGTTGAGGGGGATCAAGTTGCAGAATGATACAGATAGTACAATGACTTTTTTGTAAATTTG
AATTAAGTTTCTCTATAAACTCAGTACTATCTATTGGTAATGGATGCATATATGTATGTACATGTT
TTTGAAACGGATTGGAAGGATACAGACCAAACTCTTGATAGTGGTACCTGTGAAGAGTGGAGAAGGGA
AAATAATGAGTGTGGGAGGTGGGATGGGTATTGGTTAAAGGGGACTTCAGCTTTTTATATAAACATCCA
CTTCTCTTTCAAAAGACTTCAAGTAAATATAAGAAGATACTGATTGATTCTGGATTGTGGTATAGGGAT
ATCTTGTCTTATTTTTGTACTTTTCTGTATTTTTAAATTTCTCAAAATAGGAATGGGAGTGAGGATG
GGAATGCTGTATCTGTGAAGTCATGTTATACTGGATTCATTTCCAATTAATACTAAACATTTTATAG
AAAATATTTCTAAAATTTTATCACGGATGAATGGATGGACTTGCTCTCTATGTAATGAAATTAATCT
ATTTTATTAATTTAAAGGACAGCACAA
ACGCGTAAGCGGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_001099638.2

Locus ID:

7705

MW:

58.9