

Product datasheet for **SC214858**

KCNIP4 (NM_147181) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	KCNIP4 (NM_147181) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	KCNIP4
Synonyms:	CALP; KCHIP4
ACCN:	NM_147181
Insert Size:	1506 bp



[View online »](#)

Insert Sequence: >SC214858 3'UTR clone of NM_147181
 The sequence shown below is from the reference sequence of NM_147181. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TCCATGCAGCTCTTTGAAAATGTGATTAACTTGTCAAATAGATCCTGAATCCAACAGACAAATGTGAA
CTATTCTACCACCCTTAAAGTTGGAGCTACCACTTTTAGCATAGATTGCTCAGCTTGACACTGAAGCAT
ATTATGCAACAAGCTTTGTTTTAATAAAGCAATCCCCAAAAGATTTGAGCTTTCAGTTATAAATTT
GCATCCTTTTCATAATGCCACTGAGTTCAGGGGATGGTCTAACTCATTTCATACTCTGTGAATATTCAA
AAGTAATAGAATCTGGCATATAGTTTTATTGGTTCCTTAGCCATGGGATTATTGAGGCTTTCACATATC
AGTGATTTTAAAATATCAGTGTTTTTTGCTACTCATTGTATGTATTGATCCTAGGATTTTGAATGGT
TTTCTAATATAGTGACATCTGCATTTAATTTCCAGAAATTAATTAATTTTCATGTTTGAATGCTGTAA
TTCCATTTAAATCCATTTATATACTTTAAGGAAACAAGATTACAACAATTAACAAAAACACATAGTTCC
AGTTTCTATGGCCTTCCACCTTCTGTTAGAAATTAAGTTTATCTGGCATTTTTAAACATTTAAAATT
ATTAACATTTAAAATTAGTTTATTATCAGATATCAGCATATGCCTAATAAAAATTATTTAATAAGC
ATTTAATTTTCCATAATATGTTACAGCCAAGGCCTATATAATAATTTGGATTTGTTCAATCTTTCTTA
CAGGCTGTTTTCTATTGTATCAATCATTAGTATCAATCATTAAAGTGAAGTTGAAGAAGGCATCAACA
AAACAAGGATGTTTACAGACATATGCAAAGGTCAGGATATCTATCCTCCAGTATATAGTAATGCTTAA
TAACAAGTAATCCTAACAGCATTAAAGGCCAAATCTGTCTCTTTCCCTGACTTCTTACAGCATGTT
TATTTATATTACAAGCCATTCAGGGACAAGAAAGAAACCTTGACTACCCCACTGTCTACTAAGAACA
ACAGCAAGCAAAATAGCAAGCAAAATCACTTTGAAAGCACCAGTGGTCCATTACATTGACAACACTAC
TACCAAGATTTAGTAGAAAATAAGTGCTCAACAATAATCCAGATTACAGTATGATTTAGCTCATCATA
ATTCAGATTTATTTAATCATCTTAGCCAAAAGTAAAGTTGCCACATTACTAAAGCCACACACATCG
TCCCTGTTTTGTAGAAATATCACAAAGACCAAGAGGCTACAGAAGGAGGAAATTTGCAACTGTCTTTGC
AACATAAATCAGGTATCTATTCTGGTGTAGAGATAGGATGTTGAAAGCTGCCCTGCTATCACCAGTGT
AGAAATTAAGAGTAGTACAATACATGTACTGAAATTTGCCATCACGTGTTTGTGTAACTCAATGTG
CACATTTGTATTTCAAAAAGAAAAATAAAAGCAAAATAAAATGTTTATAACTCTA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_147181.4](#)

Summary:

This gene encodes a member of the family of voltage-gated potassium (Kv) channel-interacting proteins (KCNIPs), which belong to the recoverin branch of the EF-hand superfamily. Members of the KCNIP family are small calcium binding proteins. They all have EF-hand-like domains, and differ from each other in the N-terminus. They are integral subunit components of native Kv4 channel complexes. They may regulate A-type currents, and hence neuronal excitability, in response to changes in intracellular calcium. This protein member also interacts with presenilin. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

Locus ID:

80333

MW:

58.3