

Product datasheet for **SC214870**

PBX3 (NM_006195) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PBX3 (NM_006195) Human 3' UTR Clone
Symbol:	PBX3
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_006195
Insert Size:	1484 bp



[View online »](#)

Insert Sequence:	>SC214870 3'UTR clone of NM_006195 The sequence shown below is from the reference sequence of NM_006195. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC GGAAGTGTGCACTCGGATACCTCTAATATCTCTGGCCACACTTTCCCTGAGCTACATGCCTTGATA AGTGCATTAGAGCAATAGGAGGAAAAGGAAAGCGTTTTTGTAGCCACCATCTACAGCTTTACTGTAA AACCTTGTCTTATTCGAGAACTTGGTAAATCTGTTTTTAAAGGAATCATAATCATTTGTATTATACTT AAAAACACACAATGTTAAAAAATAAAGCACTTTATCCAATTAGGCCAAGATTTAACATTGTTGACAG TCCTGTAGCTATTTATCATAATTTATTATCAATATTTTACATTAATGGTTTCACAGTTGCCAATTACT TGGCCTTAAGGGTAAAAAGTACAATATACACTAAACCTCAACCGTTAAAGCAGATGCAAAAATTCACCT CACCTAAATTGAACCTCTTGCATATTTCCATTACTGACTTGGATTGTCTTTCTTCATATCACTAATGG AGTTGGAATAAAGAGCTGTTGCTATCCCTGTTAATGATGGTTGTGTTAAGAATCTTCTCGTCACG TTTGTGTTGATGCTCTTATGTTATAATTAGATCAGAGACTGGTAGCATCGTTTCTCTCTGAAAGCA CCAGTGCCAGAGTCTGCTCGGTAAATAAATTATGGATCCAGATTGTTCTGAGAGACGAAGATACTTGC TGCTGATAGAGGTGAAAACGAGATTGATCCGTCTGGGGTTTTACGGTGTGCACTGGGTGCTGCACAGAC TTGTCAAGGTTTGCTACGTCCTCTGGGCATCTGCAAAAGGCCCTGCTCTCTGGAGTGTGTATATAGTG TAGCAAAAGAGATTTATACATCCCAACATCAAAACACAGCTTTATTACCTCATGCAACTCATACAA ACCAATAGAATTTCAACATGTTCTGTAGCTTAGAGTGCTCACTTACTACCTCTGAACAATACTCAGCT GTAGTTTGTCTCTTTCTTATCTTTTGCATCTTGAATTAACCTTTGTTCCCTTCATAAAATGTAAT GTACATTGTAATCTTTTAAAGAAAAATCAGGGTTGCACTTGCAACTTTTAAAAACCGAGTGTGGAAA CATTGGGTCTTAATTCAACACAGGATCGGTAAACTGTTGTAATACTGAGAAACATTTTGAATGTTCT TCATCTTATTACTAATCCATGCAAAAAAAAAAAAAAGCAGCGACTAATTGTGATGCATTACAGATTTT AGTATTCAGTACTGTATATTTACCCTGTGTAATGGGGCCCCCTCTCCTTTCTCTTTTTGTATTGTA TGCGATTCTGAACTGATTGAGTCATGAAAATAATTTGTGGCGGTGATTCTAATGTATTAACCGTTT CGTGTCTCTTCTAACTGGATTACCCCTGGATTGAAAAGTCTTCTCGTGGTAGTTATATGTAGTTT CAAACATGAATAAACTTTTGTCTTCATGATTTAA ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_006195.6
Summary:	Transcriptional activator that binds the sequence 5'-ATCAATCAA-3'. [UniProtKB/Swiss-Prot Function]
Locus ID:	5090
MW:	56.8