

Product datasheet for **SC214891**

C8orf55 (THEM6) (NM_016647) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	C8orf55
Synonyms:	C8orf55; DSCD75
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PSI00062)
ACCN:	NM_016647
Insert Size:	1518 bp



Insert Sequence:

>SC214891 3'UTR clone of NM_016647

The sequence shown below is from the reference sequence of NM_016647. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

CGAAGATTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAAGGACGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
GGGCTCAGTGATGTACCAAGGACCAGTACCGCCACCTTCACACCGTCTGCCCTGGCCACCATCCTGG
GCCTGGGGGCTGCCACAGATGGGCAGTCTCAGCCATACTCTGTTCCAGCTGGAGTAGCCTCCTGACCA
GCCTGGCCACCTTGCTCCACCCACTGGGCCCCCCAGTTATTGATACCCCTCTGTGCTGGGCTCCACC
CTAGGCAGAAGGAGGAGTGGCATTGGCATCCTGACCCAGCTCTGCCCTCAAGTGGGGATGGATGGGCA
AAGGAGAGTCTGCCTGGCCCTACGATGAGGCCACATCTGTGGCCTAGGTAGGGGAGGATGGTGCCTG
GAGCAGAGGAGGACCCAAAGTGCCTCCCGAGCTAGATCCTGGCTCGGACCACTGCAAGGGCCGAGGAGCA
GGCCAGACCAGAGCATCCTGGGTACAGGCCTGGGCTCTCCAGGGCCTGGGCTGATTAGGTGCAAGTGG
GCACTCCTGAAGGTCAGAGCGGCATCTGCCAGGCAGCCCTCTGGCTTCCGCTGAGGTGGTTGCAAGGC
CTGGGGCAGAGCCTGGGTGGTCAGAGGCCGGGCTAGAGGCAGATGGAAGGGAGGCATTTGCTGACAGA
GGACGGGGCACCCGGGCTCCCCACTGCAGTCGGCCTTGCTCCTCCTCCTCTACCTCCAGTCAGGC
TGGACGGGAGGGTAGCCTTGTCCTGAGAGGGGTGAGACTAGGTGGCAGAGGGGCTCTGGAAAGACAG
CAGGCTTCTGCTGGGCGTTCCTTGTGGAGGGAATAGAGTGGGGGTGGGACTCTGCAGGGGTGTCTT
TGCCACTCGCACCCCTCGCCGCCCCACAGGGCCATGCTCTGTGACTTGGGCTGATCCCCACCTTCT
GGGCTACAGCACACAGGCCGCTGTACCCCTTAGAGCTGCCCTCTCTGGCCTGGCCGCGAGGCTCT
TTCTTAATCCTCTGTCTCTATATTACGATGTTCTTGCTGACGTGCTGGGCGCGCCTGCTTGGCG
TAGCAGAGCCTCTCTGGCAGCTTCTCAGGTCTCCCTAATGGAGACACCAGGCTACTAGGACACTGGCT
GGGGCCACCCCTCCTGCCTAATGCCTCACCTTACAGTGGGGAACTGAGGCTGGAATGGCCAGAG
TCACCAAGGCAAGTGGGGCTGGTCCAGCCTGAGGCTCCAGCTGATGCCTCAGTCCCAGAGAGGG
GGTGGCCCATCTAGCTGGGTGACAGGGGTCACTGCTTGTGACGCTCAGGGCCCTGTGCCCGCTTGCTGTT
CCCCTACATCTGTGCCTGCACATCCAGAACTGCCTCCTTGCCGCTGCTCCAGGAAGCCACCTTGAGC
CAGAGTCAAGGCTGCAGCACTGCCGATAGAACACGCCCGCCCTCACTGCTGTTCTTGCTTACAGCC
ACCATGGGAAGGCTCGAACCTTTCTGTTTTATTTAAAGAAAGCCCAACATTAAGGGTTTTTCATTGCAA
ACGCGTGAGCGCGCGGCATCTAGATTGCAAGAAATACCCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTTCGATTCCACCGCGCCTTCTATGAAAGG

Restriction Sites:

Sqfl-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_016647.3

Locus ID:

51337

MW:

53.8