

Product datasheet for **SC214903**

ZNF74 (NM_003426) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF74 (NM_003426) Human 3' UTR Clone
Symbol:	ZNF74
Synonyms:	COS52; hZNF7; ZFP520; ZNF520
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_003426
Insert Size:	1517 bp



[View online »](#)

Insert Sequence: >SC214903 3'UTR clone of NM_003426
The sequence shown below is from the reference sequence of NM_003426. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TTCTCCCTGGGGAGCAAACCTCGAAACTAACATGATGTGCTTTGGTGTCTAGTCTGCTTTCTGAGCTA
CTCAACAAGGAAAGCACCTGGTCTCCTGGCTCCTAGATCCAGACCACCTTCCTCCAGGTGTGGGAG
CCTTGCTTATCACCCCATCAGGTCTGCATGCCAGGGTGCCTCCTCTAGTTAAAGTCAGTCACCTCCC
CAGAAGGGCCACACTCCAGGAGGAGTGTGAGAGTCATTTGAGGTAGTCTTGCCACCTGTTTTCTTGA
TGGGCTGGAAGTTGTTGACAAGGGGAAAGATCTTTCTTGCCAATAAAAAAGAGGATATCGTTGGGTG
CCATGGCTCACACCTGTAATCTTAACACTGTGGGAGGCCAAGGCAAGGGGATCACTTGAGCCAGGAGT
CTAGGACCAGCTGGACAACATGGTGAGACCTCGTCTCTACAAAAATGCAAAAAATTAGCCAGATGTGG
CGGCATGTCCCTGTGGTCCCAGCTACTCAGGAGGCTGAGGTGGGAGGATCACTTGAGCCTAGGAGGTCA
AGGCTGCAGTGAGCCATGATTCACAGCACTGCACTGCAGCCTGGGTGACAAAGCAAGACCTGTGTGAA
ATTAAGGAGGTATATCAACTGTTGTATCCTCGACGGGCTCCTGACACGGCTTTAATCAGGAGTTTC
CTCCATAAACTATTATTTTCTAGATAATAATAAAACAAGAATTATGACTAGCATCACTTTCCAGTGAAC
ATTATTATATTGCTAGAGAGGAGAATAATCTTGGTGGTGGGCATTTGGAAAAAGTGAATTTCTGGAC
TTAACTCATGTAATAGTCTACTGCAGAGCTGTGTGTTAGTGACAGTGCAGTCAGGGGATCCCCAC
AGCTGTACAGCACGGCCAGCATCATTGTAGCCAGATCCTAACATGCCAACATCACCTCTTGCCATTT
AGCCCCAGTGAGAAATTGGGAGCTACCAGGCAGGTGCCAGTGCATTGAGGGAAGATGGGCACAACAT
CAGGATGGGTGTGTCTGGAAGCTCCTCTTCACTGACCAGGCTGGGCAGGGCCACCCTGGGCTGGTGAG
GCTGCCCTGCAGGGCTTCTCACTATGTAGCACCAGTCACCAGCCAGGTACAGAGAAGAGCCCTCCCA
GCACCCACTGGAGAGGGGAAGCTGAAGCCAGGGAGACGGGCTCCATGGTGGTCCCATGAGCTGTT
TTGCAACCTGGAAAAGGCGGCTCTCCCTGTCCCAACTCTTCAGAGACAGGAAGACAGAGTTAATC
TTGAATGAATGTTATTTCTACTTAGTCCTAAGCAAGAATAAGCTGCCCTCTTGGTGATCATACTTTATA
CGGAGTGCTATCGTTTATGAAATGCTTTTCACTATGCTTTTAAATTATGCAGGCAAAATAACCACTT
CAGAAAATGTGGACTTCTGAAAAATAAAATCCCCATAATTCTTAATAAACTGCCTTTTGCACCTTGA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_003426.4

Summary: May play a role in RNA metabolism.[UniProtKB/Swiss-Prot Function]

Locus ID: 7625

MW: 56.7