

Product datasheet for **SC214922**

PPP2R4 (PTPA) (NM_178000) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	PPP2R4 (PTPA) (NM_178000) Human 3' UTR Clone
Symbol:	PPP2R4
Synonyms:	PP2A; PPP2R4; PR53
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_178000
Insert Size:	1512 bp



[View online »](#)

Insert Sequence: >SC214922 3'UTR clone of NM_178000
 The sequence shown below is from the reference sequence of NM_178000. The complete sequence of this clone may contain minor differences, such as SNPs.
 Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CTGCCCATCCATCCTGTACGTCGGGCAGGAGGGGCAAGCCGAAGAGCCACCCAGGCCACAGTTCTCT
GTGCTGCTTCCCCACCCAGCAGTGGCCCTCCCATCCCCTCCCTCTGTTCTGTCCTGTTGATGAG
AGGCTGTTTACTGGGGTGGGGTGGCGAGATGGGCTTGAGGGGGCTCAGAGCATAAGGCTTCAGGGCCCA
AGTTGGGAGAAGTGACCAAAGTGTAGCCAGTTTCTGAGTTCCCGTGTGCTAGACTGGCCAGAAGAGAG
GGTCTGGGGCTGGTCACTCGGCCACTCTCTCTGTTTCTGGCCTCTTCTCCCTCACTCCCGTCCAGT
CTGTTTTTCTGAGAGCAGGGGCTGTTCTGCAGCACCAGGGAAGGAGGAGAGATACCTGTGCTTCCAT
TGCTTTTCCCTTCTGGAGTCGATGCCTTTCTAAGGGTTGGAGCTGCTCCTTGAGGGGGGGGTGAGTT
TCCAGGCCATGCCGGGTGGCCATCTATGGTAGGGCTGGAAGCTGAGGCTGGCCGCCAGCTGTGGGCT
GGGGTGGGGTGGGTGGGGTGGGTGGAGAGGCCTAGCTGTCTGGCTGGTGGCCCTCCAGGCTC
CTTTTACCCCTGCCCCCTGGGCCTGAGGCCCTGTGTCCAAGCCTCCCTGGCTCTTCAAGTTCTCTA
GCCCTTGGCTCTGCTGGGTTTCTGACTGTAGCCACATCTCTCCCGCTCCCTAAGGGTAACCTAGCCAA
TGGAAGCTGCCCTTTGGGTAGGTGCTGGGCTCCTGGGAGGGCCAGATGATGGGGTGAAGCATGTCTTT
CCAGAAGTTTCCCTGGCAGGGAGGGGATGGCAGAACTCAGGGAGGGGCTTGGGGCCATTGTATCTGG
AGAGCCTGGATTCTCTTGGCAGTCTTAGGCCCGGCCACTTCTGCTACCTTTGCGCTGCTGTGAGCCTC
ACCCTGGGCCCCCTGGGCCCTGCTTCTGCTCCCTGGGTGATGGGTGGGCCAGAAAGGTGGCAGTCCC
ACACCTTGTCTCCACCTCCCTGAAGTGTCCATTGCTTTATAGGGTGAGGTAAGAGACAGCCTCCCA
AGCCCAGGCTTTGGCACTCAGAATGGGCCAGTGGGGGCTGGGCAGGCCATTGAGGGCCACCGCCGAG
GTTTCTCTAGGGCTGTTCTGGCCTGGCTCTTACAGGCTCGTCCCCAGGCCTGCCCTTCTCCACTG
CCCCCTCTGTGCTGGGTCCACACACCCTTCAGGAAGGGGGAGCACTGAGAAGCACAGCACAGGGGCT
CAGCCTGGGATCCGGTGATGGTCTGGGCAGAGGCTGGGTGAGGAGTCCCAAAGGTGAGTACAGTTTCT
CAGAAGAGGGCCAGCGTCCACCTCTCTCCAGGGCCAGACAGCCCTTCTGGCTCCCCCATCCCCCTAT
GGGCTCCAGCCCCTTGACCCCTATTGCTGTTTCAAGTTAAAGCCTCTGTTTGCACCTGTCA
ACGCGTAAGCGGCCGCGCATCTAGATTCAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_178000.3](#)

Summary:	Protein phosphatase 2A is one of the four major Ser/Thr phosphatases and is implicated in the negative control of cell growth and division. Protein phosphatase 2A holoenzymes are heterotrimeric proteins composed of a structural subunit A, a catalytic subunit C, and a regulatory subunit B. The regulatory subunit is encoded by a diverse set of genes that have been grouped into the B/PR55, B'/PR61, and B''/PR72 families. These different regulatory subunits confer distinct enzymatic specificities and intracellular localizations to the holoenzyme. The product of this gene belongs to the B' family. This gene encodes a specific phosphotyrosyl phosphatase activator of the dimeric form of protein phosphatase 2A. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jul 2008]
Locus ID:	5524
MW:	52.9