

Product datasheet for **SC214931**

CLIP3 (NM_015526) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CLIP3 (NM_015526) Human 3' UTR Clone
Symbol:	CLIP3
Synonyms:	CLIPR-59; CLIPR59; RSNL1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_015526
Insert Size:	1528 bp



[View online »](#)

Insert Sequence: >SC214931 3'UTR clone of NM_015526
The sequence shown below is from the reference sequence of NM_015526. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
TGGATGCTGAGGGCGGAGATGCAGTCTTAGAGGCCCTGGACACCTGACAAAGAGACAGAGTCCCCACTA
GCATCTCTGACACCCGAGGAGCCCTGAGTCACCCTGAGATAGAGATTCCAGTAACACATCCAGAGTA
GAGACCCCTGTTAGCCAGCCCTCGATCATTGAGGCCCCATTATTAACAGATACTCCATAATAACCCCC
AAATACAGACCCCATGTACCCAGAAAGAGATTCCCTGAGTAGCACCTTCAGGCTAGTCCCTATCCCCA
ACCCCTCAGAGCAGATCCAGATTAACAGATTTCATATCACCCCAAATGATGGTGACCCTCTCCACA
TAATGCATTACAACAGAACATTCTGAATCACCCAACCCTGGATCAGAAACCTCCCCATTAACAAACAC
TGCCCCCTTAAGTCTCTTGAAATAAACATAGGTACACCCCCAAAGCAAAAGAGTAACAGACATTATG
TCATTGTTCCCCATTTAACATCAGTCTCTCAAGATGTCGTGACCCCATGGTCACCCTGAAGCCCTTAG
ATTCCAACCCCTCAATCAGAGACTTCCTTCATTACAAAGACCCTTGTTCTTATCCCTCAAGAAGAAAC
CCACCATAACCAGCCCACTGTACCCCTAATTTACAGACACCAAAACAGTCTGGAAGTGCTAATTACA
GGACCCCCCAAGTCTTCTACCTCTGCACCCTCAAGAAACCCCAAGTGCCTTGTATGAAGCCACCCC
ACATGGCCACAGCTCCTGTGCTGGCCAGACTCCAGAAAATTCTCTATTTTTTAAGTAACGACTTCCC
CCTTTGGGGGACCCCAAAATTTGGAGGCCCATTTCTAGGACTCTGGGGATCCCAACCCTAGAGTACAC
ACGTCCCAAACTCCCCTGTGCCCTCAAGTCTACAGCCCTAGAAGACCCCAATGCCGTAACCTCTAGG
ACCCCAAAATCATGGAATCCCAATCCCAGGGAATCCCAATTTGAAAATCCAATCCCAAGTCCCCAG
GAAACCAATCATGAGGTCTTTGTCCTGGTATGGAGGAGACTGCAGTCAGGATATGCATTCCAGGCTC
CCAGACACCTCAAGCCCTATTACAGGCACCAAGGAAACCCACACAGGAATTCCCATCCCTGGAACCTG
GAGAATTTCAATGCCCGAGTCCATGGGTTTCAAGACACCAAAATTCAGAGCCCGAGCCCTAAGGGAA
CCCCAAATCCTAAAGCCTCCATCTCTAATAAATGGAAGGCCCAAGGCCCTGAGGGGATCTCAAATCCT
GGAACCCCGATTCAATCTACGTTCTAGTCACTGGCCTCAAAGGACCCACAGCACCTGGGCCAGACCA
ACAGCTCGAGGGAGAACCTGAAGGCCAGGGGTCCAGGGCGGACCTGGGGCCCCGACCACCAAGGACA
GCTCAGACTGCCCTTCACTGCATGTCCCAAACTCAGCATGACTCCTGTCTTCAATAAAGACGT
TTCTATGGCC
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_015526.3

Summary:	This gene encodes a member of the cytoplasmic linker protein 170 family. Members of this protein family contain a cytoskeleton-associated protein glycine-rich domain and mediate the interaction of microtubules with cellular organelles. The encoded protein plays a role in T cell apoptosis by facilitating the association of tubulin and the lipid raft ganglioside GD3. The encoded protein also functions as a scaffold protein mediating membrane localization of phosphorylated protein kinase B. Alternatively spliced transcript variants have been observed for this gene. [provided by RefSeq, Dec 2010]
Locus ID:	25999
MW:	56.7