

Product datasheet for **SC214934**

FLJ44635 (NM_207422) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FLJ44635 (NM_207422) Human 3' UTR Clone
Symbol:	FLJ44635
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_207422
Insert Size:	1526 bp



[View online »](#)

Insert Sequence: >SC214934 3'UTR clone of NM_207422
The sequence shown below is from the reference sequence of NM_207422. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GCTAATTTCAAAAACCTACCAAGAAACATGAATCCAGATGGCATGGTTGCTCTGCGGACTACTGTGAGGA
TGGTGTGACCCGATATATGATTTTCTTTAAGGATGGTTTAGAAATGGAAGAAATGTTAACAAATTTGGCA
ATTACTTTGGATCTATCACCTGTCATCATAACTGGCTTCTGTTTGTATCCACATAACACCAGGACTTA
AGACAAATGGGACTGATATCATCTTGAGTTCTTCATTATTTTACTGATTATTTGGAGTGGAGGCATT
GTTTTTTAGAAAAACATGTCATGTAGGTTGTCTAAAGTAAATGCATTTAACTTAAAAAAAATACC
TGAGACTGGGTAATTTAGAAAAACAATAGATTTCTTTATCTTGTGGTTCTATAGGCTGAGAAGTTCAAG
GGCATGGCCCTGGTTTCTGGCAAGGGCTTTTGTCTGTGTCAACATAGTGGAGAAGATCAAAGAGAA
AGTGGACACATGCAAAGAGGGGAAAACCTGAGTGGCATCCCGGCTTTATAACAACCCACTGTCACAGGA
ACTAATCCATTCTGAGAGAACTAATCCAGTCTTGTGAGACTGAGAATCACTCACTCCTGCAAGACTG
CCACCAAGCCATTACGAAGGATCCATCCCATGACCCAGACACCTCCCACTAGGACCCACCTTCCAAC
ATTTGGGGATCAAATTTCAACATGAGTTTTGATGGAGACAAATAAACACATCCAAACCATAGCACAGG
GACTGCTAAATACCCCAACAACAGAAATAGAAAGCTCCCTAAACACCTTCTTTAACTGAGACCTGGCT
GTCTCCTGAGAACCCACTTCGCTGTAGGCTCTCCTCAGAGATTCTTTCCCTCCCGCAAATACTTCT
GCATCTGGAGGTGGAGTGGGCCCTTCTTGTCTCCCACTGCACCTTCCAGACCATGACCACTGTACCT
CTCAGGCACTTCCCTCCATTCTATCCTTTTGGCCACCCCTAGTCTTGCCCAATCCAAATCACTTCAA
CATTGATGTGAAAGATCTATCCCACTTCCCACTTCTTGTCTTCTGGCCTACTCACCTCTGTTGCCC
TCCCTTTCTGCTCCAGTACCCACTCTCAGGGCCCAAGTGGATCTTGCGATGGCCGAGAACTGCTGCA
CCATTGAAATCTTAACCTCAATCCCAACCTCACATCTTTCTAACTTTCCCATCTGATTACTTCATT
ACACTGTTGAGCTTTGTTGAGACTTGTATGAGCTCATCAGAAAGAAATCCCAATGTCCAATAAAC
ATGGAATGTGCTTAACCTCATCTACTGTCAGGAAACGAAGAATGAGGGAGAATGGGGCATGAGGAAG
GAGCATGAGAAACTTTTGGAGTAATAGATATGTGCATTATCTTGATTGTGGCTCTGGTTTCAAGGAT
GTTTACTAATCAAATCGTCACTTAAACGTGCAGTTTGTGTACATCAATTATATCAATAAAGTTGT
AAAATATA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_207422.3](#)

Locus ID: 392490

MW: 57.4