

## Product datasheet for **SC214949**

### **FAM13C1 (FAM13C) (NM\_001001971) Human 3' UTR Clone**

#### Product data:

|                |                       |
|----------------|-----------------------|
| Product Type:  | 3' UTR Clones         |
| Symbol:        | FAM13C1               |
| Synonyms:      | FAM13C1               |
| Mammalian Cell | Neomycin              |
| Selection:     |                       |
| Vector:        | pMirTarget (PS100062) |
| ACCN:          | NM_001001971          |
| Insert Size:   | 1511 bp               |



**Insert Sequence:**

>SC214949 3'UTR clone of NM\_001001971

The sequence shown below is from the reference sequence of NM\_001001971. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGCAAGCAAGATGTGGCCAAACTATTGAGGTTCAAGAAATGTTATGATCACTTTCACCCATGATATA
AAGTAAAGTTTATTTCTCTGCCATCCTTGCTAAGTAGTTTTGACACAATGAAATGGAAGCACTTTA
GTGGTAGTATTAGCTGTTTTAAGAAGGAATAGCAAGTTTAATTATATACAAGGAGAAGGGATTTAAAC
GGGGGAAGAATACAACAGGTAGCCATATAATTGGGAAAAAATTCAGTGCCTCCATGCCAAGCAGAAA
ACTCATAGTCAATACAAGTATTTTTAAAAATGTCTAATTTTTATCAAACTAAATAACATAGCTAGGA
CACTTGTTAGGAAAGTTTATTTAGTATCCAAAGACTGTTTATGTTGATGTATGAAAAGAGCATGATT
TAAAAAATCAATCATAGGAGGAAAAGAAATTCGCTTTTCAAGTAGGAAGGAATACAGCTAGCAAGAAA
GCAATTTATTTGAACTTCTAATGGATTTTGTAGTGATAAACATTTACTACCTTGCCTTTAAGTCTG
CTAGGCTCTCAGTACCCTAAAAATAAAGTAGATTGTGTTGCTATTTTTTTCTTTCTCTATAAAAAATAACA
CATTATTTTATCCGTTATTTGAAATTTTACATTTCTGGTTACCAAAGTTCATTCTGATAGCATGTAATT
TGTGAATTATTATCTTTGTCTATAACTGACAGATGTTTATATTAATAAATAATTTGTATTAATAAATTT
AAAATAGGTATTTGGATAGATATGTGCTGTAGTATATAATCTAATGTGTCCATAGTATTATTGCTAA
TCTTTTGGTTTACTATAAGATGATATAACTATTTTTTTCATTGGGAATATACATTTTCTTAATGTTCCA
ACATCTATACTTTGTAAAGTCAAAACATTTCCCATGAGCTGTAGTTATTCATCCTTCTGTACAAATGA
AAAGTTTGGAAATGTTTGGCCTGATACCTTGAAAAAGAAGCCAGAATATTTATTTGCTTCATCAACTT
CAGTGTATATCATTTTGTGTTATTTTATACGAAAACATGTTTATTATTTTCATTTTGTAAAAGGAAGT
AAAAGGTCAACATTTTCTCTCATGTACCAACCTTGTTTGTATTTCTATTTTCTGTAATGTTTAAGTAGT
ATGTTGAAGAAATTCACATTCTCTTATAGTTTGGATGGGAAGACTATTGACTATTTCAAGAACAGACTT
ATTTACAGAGGCTTATTGTTTCTCTGTATTTACCTAATATTTTATAACTTTTATGAATCAGAATAATGT
CCTTCATAAATTTGTTTAATTGAAGTCATCTACTTCTAACAGGACAGATACACAATTTGAGGTTTA
CAAATTACATCTTTGATAAGGGAAATGGTTTCGTGACATGTACACAGTTGCTATTAATAATGTAACCTCA
TATATTCTATATGATTGTAAATATTTTATACAACAATACAAATAAAATATTTTCTATTATA
ACGCGGTAAAGCGCCCGGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

**Restriction Sites:**

SgfI-MluI

**OTI Disclaimer:**

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

**Components:**

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

**Note:**

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

**RefSeq:**

[NM\\_001001971.3](#)

**Locus ID:**

220965

**MW:**

58.9