

Product datasheet for **SC214988**

NCKAP5 (NM_207363) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	NCKAP5
Synonyms:	ERIH1; ERIH2; NAP5
Mammalian Cell	Neomycin
Selection:	
Vector:	pMirTarget (PSI00062)
ACCN:	NM_207363
Insert Size:	1531 bp



Insert Sequence:

>SC214988 3'UTR clone of NM_207363

The sequence shown below is from the reference sequence of NM_207363. The complete sequence of this clone may contain minor differences, such as SNPs.

Blue=Stop Codon Red=Cloning site

```

GGCAAGTTGGACGCCCCGAAGATCCGCGAGATTCTCATTAAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC
AGCGCTGCCCCAGAAATTGAGACAACCTGAGAAACAAAAGACGATCCCAGAAATAGATTATCGAAAAT
TTCCCTAGAGTCAATTAATAATTTAACAGCAATACTGTGATTTTATTAGAAAAAGAGAAGAACTCTCT
GAACAAGGTTGAAGGACAGAAGGAAGAAAAAGAAAAAATGAAGAGACATCTTTGAGTAGTTCAGATAG
GCCTGGGGTAGACAACCTGGAATCTTTGAGTGATTCTTTATATGATAGCTTCTCTTCTGTGCCAGTCA
AGGTTCAAATGATGTATAAAGGACATCTTCCCTTAGTGAGCTGGGACTGGAGCGCTTAAGAAATGAT
GGGTGGGGGGTGGGGGTGCACCGCTTGATAGAGATAACAATAAACTATTGCAGTACCAGAGCCTTCTCT
TGTCAAATTCACAGCAGGCAACCCACCAGAGCTTATTTCTCTGACAGGGCAATAAGATAGACTCCATT
TATTGTGTTTCAAGAGGATTAAGCGTAAACACATCTATGATACAGAATCCTTAATTTTGCACTTTTTT
GAATATTTGTACAGAAGTTGTAAATTTTTTGGAAAGAGAAATTATTTGTAGCAAAAAAGACAGCAAT
AAATGGAATCAGTGCCATGCTCTTGAAATAATGTAAGTCTTAGAAGTTGATGATAATATATATTTT
TTAAATCCCAACTGAAGTTTTTGTGAAGTTCTGTCTGGTCTCAAAATGTTTGTGGGTACACTCT
GTAAACCTACAACAGGGCTGCCAAAAATCGGAGGGTTCCTCCTCATCTCCATCTCACAATCTCAAT
TTGATGGAAATGTTCAATTTTAGTGTAATTCAGATTCTGTCAGAGATTCAAGTGATAGTAATAAGTG
TCATTCTGCTTCTGCTGAAAAATGAAAAGGCTCCTGAAGTGTGGACACTGATTGGGAGTGTGACATTGT
ATCAGAAATGACCGAATCTATTCCCAATACCAGTTTTCTTCCAGACATTTCTTTAGATTGTCTTTT
ACTTAGTGCTTCTCTATGATCTGAATATTATTTGATTTTTATCTTCTGCTTTTTATTTAACTCTG
GGCACTCTAAAAATGAAAACAAATTTCTATTTGCAATGTTCACTTTTAAAAATAAAAATTAATGGTGCTA
CGAAGAATCTTTTTAATATCTTTTTTTCTACAAAGACTGTTTATATGAAGGATAAATCTATTTT
AAAGGTTATGTGTATTTTTCTAGATGTGAACATTTATAATTACTTATGTACAGGAGCTTGAAACTA
GGCCCAATAGAAAATTTTTAGGATCTATATGGCTACTTTAGCACATAATTGTTTCTTTAAAGAGTATT
GTATGATCAGTGTTATTTGGTTAATTTGTGCAATTTGTTTTATTTTATCTTAAATGAAAATTATGTAAA
ATGTCCTTGTCTTTCAGACTTTAAAAATCTTTTTGTTTCTTTCTGAATAAAAGTTATATCACATTTG
CACCTCTGAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
  
```

Restriction Sites:

SgfI-MluI

OTI Disclaimer:

Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components:

The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

Note:

Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq:

NM_207363.3

Locus ID:

344148

MW:

60.6