

Product datasheet for **SC214995**

SPOP (NM_001007230) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SPOP (NM_001007230) Human 3' UTR Clone
Symbol:	SPOP
Synonyms:	BTBD32; NEDMACE; NEDMIDF; NSDVS1; NSDVS2; TEF2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001007230
Insert Size:	1539 bp



[View online »](#)

Insert Sequence:	>SC214995 3'UTR clone of NM_001007230 The sequence shown below is from the reference sequence of NM_001007230. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCCCACGCAAAACGCCTGAAGCAATCCTAAGATCCTGCTTGTGTAAAGACTCCGTTTAATTTCCAGAAG CAGCAGCCACTGTTGCTGCCACTGACCACAGGTAGACAGCGCAATCTGTGGAGCTTTTACTCTGTTGT GAGGGGAAGAGACTGCATTGTGGCCCCAGACTTTTAAACAGCACTAAATAAATTGGGGGAAACGGGGG GAGGGAAAATGAAATGAAAACCTGTTGCTGCGTCACTGTGTTCCCTTTGGCCTGGCTGAGTTTGATAC TGTGGGGATTGAGTTTAGGCGCTGGCCGAGGATATCCAGCGGTGGTACTTCGGAGACACCTGTCTGC ATCTGACTGAGCAGAACAAATCGTCAGGTGCCTGGAGCAAAAAGGAAAAAAAAAAGAAAGGACATTG AGTTTTAACAGAGGGGAAAGGAAAGAAAGAAAGATTTTTGCAGAATTTCTCAAAATCAGTTTGGA TTCCAGTAGTATTATATTGAGAGAAACAAATTTTAGTCTTCTAACTGTGCTAAACTTGGATATTG TGAAAACCTCTTACCACCATACAAGCATCAGAAGAGCTCTCTTGTGTTAGCACTTATTGTTTGCAAGA ACAGAATACATCCTTTTATCCTTTTATGAAAAATGACAAGTGAAGGCAAAAGGGGAAGGTTATTTGATC TGGAAGATGAGTGTCTGATGTGGTGGCTTTTGCAAAATCTTTATTGGTGTGAAAACCTGGAAAAAT AACTCATCCGAATTCATATTGTCTTGACAAGAACTATGGTCTCTGTTTTAGATATTGTGAAAAATG TTTTTGGGCATTTTTCTGATTTTATTTCTTCTCCCCACCCCTTTTTCTAAAAACAAACAAAAAA AAAACACACAAAAACAAAAACAGAACAAAGAAGAGAGAAGGAAATTTTATCAATTAATAATGCTGTGTG ATAAAATCCCAGCCCAGATTGCTCAGCTGTTTGTACCTGACTTGCCGCCTGCATAGGAGCCAGTTCTGT TCCTTCTGACTAGCCCTCTTCTCCAGGGGAGAACTTCAAATGTTAATTTTTTTTTTTGAAAATA TAAATAATTACTATTTTGTACTGTGTGGTATCTCTGGTCTTTTGTTCCTCACCTGCCTTGTCTCTTG GGTCTGAGTCCCTTGCTTAAGGGATTTTGAAGTCTAGTTTTTCACTTGCAGAGATTATGTCTGAAATG CCTAATGAGTCGAGGGATTTGTTGAGACTCCGTAATCTCAAGTCTCTTTGTGAGCTATCAGCATCTG CCAGTCTCTTGTCTCCCTGAGTATCTCAGTCCATATCCTGATGAGGGATCAGGCCCCACCTCTGC CAAGGCAAGTAATGGTAGTGGGCTTTTAAACTGCCCCCGTATGTTTAAAGACCTAATCCCCACCTCCC TTCTTCTAACTAAATATAAAAAGATCCAGGGGACATAAATGTGGAGATTAATAAAGGGAAATTATTGT CTCTAAAAAAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001007230.1</u>

Summary:

This gene encodes a protein that may modulate the transcriptional repression activities of death-associated protein 6 (DAXX), which interacts with histone deacetylase, core histones, and other histone-associated proteins. In mouse, the encoded protein binds to the putative leucine zipper domain of macroH2A1.2, a variant H2A histone that is enriched on inactivated X chromosomes. The BTB/POZ domain of this protein has been shown in other proteins to mediate transcriptional repression and to interact with components of histone deacetylase co-repressor complexes. Alternative splicing of this gene results in multiple transcript variants encoding the same protein. [provided by RefSeq, Jul 2008]

Locus ID:

8405

MW:

59