

Product datasheet for **SC215004**

ZNF543 (NM_213598) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF543 (NM_213598) Human 3' UTR Clone
Symbol:	ZNF543
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_213598
Insert Size:	1534 bp



[View online »](#)

Insert Sequence: >SC215004 3' UTR clone of NM_213598
The sequence shown below is from the reference sequence of NM_213598. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

CACCACTGAAGAAAAATCTGTGG**TGA**AAGGGAACATCTTACCATCTGGCCATTACACTGAAGAGAAACTT
CATAAGCATCCTCTCTTTGAGAAAAACGTGTAATAGATGATCATTTGTCGTCTAAACAATCAAGTTAGAAA
TTGAACCAGCCTGGAGTCTTATTCTCCACCTGAGAATTCACCCATGAGAGAGACCCAGTGGTTGCTATGC
ACTTAGAAAACTTTAGCCACATCTTTCTTCTAGTTTACAGTGAAATATTATCTCAGGGATATTTTTT
TTTTTTTTTTGAGACGGAGTCTCGCTCTGTGCCCCAGGCTGGAGTGCAGTGGCGGGATCTCGGCTACTG
CAAGCTCCGCTCCCGGGTTACGCCATTCTCCTGCCTCAGCCTCCCAAGTAGCTGGGACTACAGGCACC
CGCCACTACGCCTGGTAATTTTTTTGTATTTTAGTAGAGACGGGGTTTACCCTTTTAGCCGGGATGGT
CTCGATCTCCTGACCTCGTGATCCGCCCGCCTCAGCCTCCCAAAGTGCTGGGATTACAGGCGTGAGCCAC
CGCGCCCGGCCATCTCAGGGATATTTAAACAAAGGAAGAGGAAATATAGGGAAGAGAAAGAAACCACAGC
TTCTTTTAGACTTGCTGACAAGCCTATGGCAATTTTGTCTCCCTTTGTTATTTTATGTGGGAGAGGGA
AATGTTTCAGAAACGAAAGATCTCATCCCTTTATTTTCTGTGTATGCATTCAGTCTGCCAGTGTC
GAACAGTTATACAAACGAGGGCTTGAAACATTTTGTGAAACTTGCTTAGTCTGCAACGTTCTCCACT
CTCGAGGTGCAGAAGCATACATCTTTATAAAAAATATGGAGGCTTGAGTTATGAAACACTTTTATAAGGA
GATAAGGATGTACATATGTTTGGACACATTGTCCACTTCAGATGCTTGTGTTAAAAAATGTTTACCATG
TGCTTTTACCCAATACACATACTTTTCTTGATGTGGATTACTTTTTTTTTTTTTTTTAAAGTTTCTCT
GGGATAGGGATGATAGTGTAATTGTTTCATGCCACAGTAATCCTGGGAATCTTAGTTTTTCTCTATGATT
GTGGTTTCTCAGCTTCTTTACCTCTCCTTTGAGTTAGGTCCACCTCAGTTTTCTAAAGGTCAAATTAAGG
CCCCCCCCAAAAATACAGTATGTAGCTGAATTCCTTTGAAAGATCATGATAAACTCAGAAGAACAAACATA
GAAGGAATGTCATTTTAAATCTTTTAAAGAGTGTGTGAAAAATCACTTTAACTGGAATTTTAAACAGGCA
AGCCGGCATGCATATATGTAGGGATGTGTGACTTTCACAGTTGTCACCATGGACTTAGGCTTGATTAAAA
TGCTACACATTTAGAGCGTTTTGTAATTCCTTTTTATATTTTAAAGACTAGATTTTGAAGCCATGTGATA
GTTTTACCTCATATTAATATATTGTACAGAAGTAAAAATAATCATGTATTCACTTTGAAGTCC

ACGCGTAAGCGGCCGCGGCATCTAGATTCTGAAGAAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_213598.2](#)

Summary: May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]

Locus ID: 125919