

Product datasheet for **SC215006**

ANKRD26 (NM_014915) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ANKRD26 (NM_014915) Human 3' UTR Clone
Symbol:	ANKRD26
Synonyms:	bA145E8.1; THC2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_014915
Insert Size:	1504 bp



[View online »](#)

Insert Sequence: >SC215006 3'UTR clone of NM_014915
The sequence shown below is from the reference sequence of NM_014915. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
CAGGTTTTAAAGAAAAATTATATGATCTGAAGATAAGATTTTATTACTGGGCTGTTTATGTGACATTT
TGTTTCTCATTAATATTAATATATGACATGTAAAAATATTAAAGGAACGTATTTTGTATCATATATG
TCAGATCAAAATTTATATAGTATTTTAAACAATGTTACTTGGCTTCTTTAACTTTTAAAGTGGATTTTGC
CAATGAAAACAGGAATATTGAGTTTTACATACTCAAAGTCCCAATGTGAGCTGTTTAAACAGCCAA
ACCACTAAGTTGTCATTAATACTTTAGTGCAGTTAATTGATGGCTTATTTGTATTTTGAATTTTATC
ACTATATGAGTGGTAGATGTCTAATTGGACTTAGACATCATTTTGTGTGCAGGAAGTGGTCATTG
TTCATGTTTGGATGATTACAATTGTTACAGTGCCATCAAATTTACATAATTTTAAATTTTAAAGCACTGA
TTTATATAGCTGTTCTCAGGGAGAAACAAGATTTGCCTGATTCACCTTTTCTGTTTGTATTTATTTT
TAATTAACCTTTTGGAAATAGCCTTCAGTTTAGAGAAAATTTTGAAGTTTATGACAAAGACTTTCTCC
TCCGAACAATTTGAGAGTAAATTGCTGACATTATGTCATCACCCCTGATTACTTTAGTGGGCACTTCTA
CAAGTCAGGACATTCTACATACTACAGTGCAACCTTGAATCAGAAAGTCAACATTAATACATTCTA
GTTTGTACTTCTCTGTAATGTCCTTTATAGCAAAAGGATCCCATTCTAATCCATGGCTTGTATTTAGT
TGTTGTATCTGTTAGATCTTTCAGTCTGGAACAGTTTCCCTCAGTCTTTCCTTGACCTTCATGACCTTGA
TATTTTGAAGAGTAAAGTCTGGTAATATAAGAAATGTCCTCAATTGGGTTTCATGGTGATTAGATTTAG
ATTATATGATTTTGGCAGGAGTATCACAGAAGTGACCCAGTGATCTTCTCATTGCATCTTACCTGTG
ACCCTACAGTTTGTATTTGCCCTGTTACTGGTAATGTTAGCATTGATCACTAACAGTGGTATCTGCTGA
ATTTCTCCTCTGTAATGTTATTTTTCACCTTGTGGTTAATTAGTATTTTACAGGGAGGTAATTTGAGA
CTATGTACATATCTACTCCCCATCGAACCTCAATTTATTCATTTGCTTGTATCAAAATAAACCATGG
ATTCATGGTTTATTTTATTCAGTGAGTTATAATTCATTAATATTGTTATTTAAAAATTTGTTTCTTTT
AAGAAATGTTATTTTGAACCTTCAGACTTACAGAAAGATTATAAGAGGTAGTACAAATAATTTTGTG
GATATTATTTCCCAATGTTAACATTTTACTGCATTTACTTTATCCTTTCTCCCCCTTCTCCTCTGTC
TTTCTAGATGAATATGAATATAGGTAATACAGATTTTTTTTCTAACTGTT
ACGCGTAAGCGGCCGCGCATCTAGATTGGAAGAAATGACCGACCAAGCGACGCCCAACCTGCCATCA
CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_014915.3](#)

Summary: This gene encodes a protein containing N-terminal ankyrin repeats which function in protein-protein interactions. Mutations in this gene are associated with autosomal dominant thrombocytopenia-2. Pseudogenes of this gene are found on chromosome 7, 10, 13 and 16. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2011]

Locus ID: 22852

MW: 58.7