

Product datasheet for **SC215007**

CBLN2 (NM_182511) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CBLN2 (NM_182511) Human 3' UTR Clone
Symbol:	CBLN2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_182511
Insert Size:	1518 bp



[View online »](#)

Insert Sequence: >SC215007 3' UTR clone of NM_182511
The sequence shown below is from the reference sequence of NM_182511. The complete sequence of this clone may contain minor differences, such as SNPs. **Red**=Cloning site
Blue=Stop Codon

CAATTGGCAGAGCTCAGAATTCAA**GCGATCGC**

GGCTTCTTGGTGTTCCTCTA**TAA**ACACAGAGCCCCCTAGATGGTGGGGGAATGGCAAACCTGGACCCAGG
ACTCCGCCCTTTAAACACCCCTGAACCTACTGGAATTGGACACCTGTTTCCAACCTCCGTCAGACTGTT
GCAGTAGAAGAATGATTTCTTTGAAACCTCCAGTACTTTGTTTTGTTTTTGAATACTGACAATTC
CTCGGGAACCTGGCCTAATTAGTTTTAGATGACAAGGTCTTAAGGAGAAATGAAATTATCGATTTGAG
CAATTTGTACCTGTGATTGTAAAGTCAATATCGGATTTTATTGTTGGGACCATGGACCTCTTTGTTTGT
ATGTTGTATTGTCGTCCTCAACGGAAGGAGAGCTCCTGACTCCAGGATGGGCTGCAGGTTGCAGTCAGGGC
TTGAAGTAGGAGCCCAGCAAAGAACCACCTGCTGGACAGTCCTTGACATGTGTTCTGTGTGTCTGTAT
AGCCTTAAGAAAAAGAAATGGCTTCACTTTCATTCTGTATTCTTCCCCCACCATGTGGCTGGGAGGACTT
GGGAGGGGGATGGGGACATTGGGAACCTGTCAAGAAAGTCTTTATCCAGAGAAGCAAATTTTGCACGATT
GGACTGCAATTTTGTGTTGATTGTTTGTGTTTTCTTGAAAAGCTTTACTTTTCTTCCACACTCAG
CTCTCCCTCCTCAACCCCACTTTTATTTTCTTGCTGGGGTTGAGGAGAGAAAATATAGAATTCCTGGAT
AAGACCAAAACAAAACAAACATTAAATACCTGTATGTTTTGTTTAGACGAGACCAAACTAAACAAAA
GTATCTGTTTATCAAAGTAAAGTAACACAATGGACAATTCTGCTTATTCTCTCAAAGAGATTCTAAGAT
GCACCTTTAGAACTATTAATAGCAACCTGCATTTTTTTTAATTTATACTTCAGAATCCTTTAAGAACCT
GGTGTTCCTGAGTGGTCTGAATCATATAAGTTGGTAATGGAAGCTGTAATGACCAAGTCCCCTAAACAT
ACTATGTCTTTGCCACGTGTGCTGTGACTTCTCTGTGGGTGATTTAATTTATTTGGATCCACCTCTGAGT
GAGCGCACAGTGATCAGGTGCTTCAAAGCCAACAGACCAGCTCCTCTCTCCGATCCTCTTTGATCT
GCCCAGGAAAGGGATGCATTGACACTCTCCTGCATGCACCTGGCGAGAAGCCACCTGAAAGTCACTGTGG
TTAAAGATATTGGTGGAGGTACCCAGGAGCACTGTTACAAATCCTTCTGTTTTGGCATCTCGTACAAC
ATTATTAAGACACAGCTGAGAGTTGATGGGTGTGTAATGCATATGCCAAGGAAATGCTACTAATCCCAA
GCAATCAAAAAGGAGACCTCAAACCAGATGTTAATTTGTTCTTTGTGTAACAATGAACCAAAATATTGA
TGATAAAGTCATAATTTAAGATTCAGAATAAATGGGTTTGATGTCTG

ACGCGTAAGCGGCCGCGGCATCTAGATTCGAAGAAAATGACCG

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_182511.3](#)

Summary: May play role in synaptogenesis induction.[UniProtKB/Swiss-Prot Function]

Locus ID: 147381