

Product datasheet for **SC215014**

C5orf56 (NM_001013717) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Symbol:	C5orf56
Synonyms:	FLJ21568; FLJ26121
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001013717
Insert Size:	1543 bp



Insert Sequence:	>SC215014 3'UTR clone of NM_001013717 The sequence shown below is from the reference sequence of NM_001013717. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GCTTATAATCAAGTAGAGAGACAGACATAAATGATTAATAGATGCTCCTGTAATGGTGGTTGTGAGGTG TGCTGCCACAGCCAGCCACACACGCCAGTTCCCTTCTCTGCACATGTACATGTCTGCTCATCCTTCCC ACCTGAACTCATGTGCCACCTCACTTATGCCCCCTTCCCAATGCCGACCCACACAACAGGGTGTCCCTC TGTGCCACAGCAGCCCTTGGGTATAACTGTTCCCCACAGAGGTGAACCCCAAGGGCAGAGTCCAGGC CACACTCACCACGGAGAGCCCTGCGCTGGCACAGAGCCAGCCAGCAGAGATGCTCAGAAGATGAATGA AAGTCCTGCTGAGATTGATCCAATCTGTTTGTGAACATCTCTGGGGTCCAAGAATTCACGGCTTCTT TTTCTTTAGAGAGCCTTTTGTTCCTTTTGGCATTATATACTATATTTGGGTTAGGGCACAACAATGTG AATATACTTAATACTGAACTCTACACTTGAAAAATGGTTAAGATGGTAAAGTTTTAAGTTATATGATTT TGCCACAATTTTTTAAATGCAGAGCATAGGTACAAAGAAGGAAAAAAGGGGCAAGGTCCATTTTAT TGAGAAATAAAGAACAAGGGTTCAATTCAGAGAGGATTGGATTTTAAACGTAATTTAGATTAGGGAA ATGCAAATGAAACCATATAAGTAAGATACCACTTCACATCTGTTAGGATGGCTATTTTAAACAACA AAAACAACCTGAAATAAGTGTGGCAAGGATGTGAAGAAATAGATACCTTGTGCATTGCTATTGGGA ACGTAAAAATCATGCAGCCTTTTGTGAAAACAATGTGGTGGTTCCTCAAAAACCTTAGTTATCATATGATC CCGCAGTTCCACTTCTAGGAATATACCCAAAAAGATTGAAAGCAGGGACTCAAAGTGATATTTGTACAC CAGTGTGTGTAGCAGTATTATTCACAATAGCCAAAAGGTGGAACAACCCAGATGTCCATGACAAGAGG AGTGGATAAACAATAATGGTGAATATGTACAATGGAACATAATTACGCTTTAAAAAGGAATGAAATTC TGGTACATGTTATAACATAGATGAACCTTAAAAACATTATGCTAAGTGAAATAAGCCAGACACAGAAGG AAAAATACCAGATGATCCCACTTATATGAGGTACCTAGAATAGTGAAATTCATAGAGACAGAAAGTAGA AGGATGGTTGCCAGGGGCTGGAAGGAGGAAGAACAGGGAGATACTGCTTAGTGGATGCAGAGTTTATGT TGAGGGTGATGAAAAAGTTCTGGAGATGGATGGTGGTGGTGCACACAGTGAGAATGTACTTCAT ACTATGGAACTCTGCACCTTAAAAATGGTTAAAAATGGATGTAGACATTTTATGCTATTTTGGCATAACA AGAAAGTAATTTGGGTGTATTTAGAGATTTTAGAAAACTGTAAAGCACAAGAGGAAAAAATTAATA TCATTGCAATCACCAGGTAGAGAC ACGCGTAAGCGGCCGCGGCATCTAGATTGAAGAAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	NM_001013717.2
Locus ID:	441108
MW:	60.1