

Product datasheet for **SC215035**

SKA3 (NM_001166017) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	SKA3 (NM_001166017) Human 3' UTR Clone
Symbol:	SKA3
Synonyms:	C13orf3; RAMA1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001166017
Insert Size:	1538 bp



[View online »](#)

Insert Sequence:	>SC215035 3'UTR clone of NM_001166017 The sequence shown below is from the reference sequence of NM_001166017. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CCAACACAGAACTGAACAAAATGAGATGAAGCCGAGCTGGACCGATTTTAACATTACATTGCCCTG CCTCTGTCCCCCTTTAAACGTTGACCCATTTTAAAGACAAACATGAACATTAAACATCATAATATGCTTT TTATGAAGTTTCAATAAGGTTTAACTTAGTCTTGTGACATGTAGCCAGTCATTCACTCTTTAAGGA TTATTAGTGTTCATTGATACTAAATTACCCAGCTTAATCAACAGAATGGTTAAGTAGTACCAGGAAG TAGGACAAGTAATTTCAAAAATATAAAGGTGTTTGTCTACTCAGATGAGGCCGCCCTGACCTTCTGGCC AGAGAGACATTGCTGCCAGCCAGCTCTGCCTTCCCATCATCTCCTTTAGGACCGTCCCACACCTTTTA CTTGCTCAGTGCTGTGAAGATGCAGTTGCTGTTTGCAAACAACAGGAACACCAGTTAACTAATTAG GAAAAGAGGGAGATTTCCAGGCCTGGGTAACATATACTGTGACCATTTGAGGTAGAGACAGGTCTCAA CAGTTGGAACCAGGAACCTCTGCTGTGAGTTGAGAGTTTTGTTTCTCTCCAGCTTTTCACTGTGTGGG GGTCTTTTCTCTTATGTGAGCTCTTCTATCACATGGCAGCTGACCTCTCACGCTCCACTCTGCAGCTT GGACACCCAGTAGACCTGAATTTCACTCTCTCTAAAAGGTTCTGAGGGCTCATCCTGGGCCAGGGGCC CTCCTGTGCACTGTAGCTATGGCCACGGGAGCCTCCAGAGCTGCCTGGTAGCTTCAGGTTGACCTGCT TATCAGGCCTACGATCCTTCTGATTTAAGTACAGCTGGAAAGTATTATCTAATTAAGTTCATGATAGTG CTTTTGGAGAAGTTGTCAAATTACAGCCAATGAGAAAATAAGGACCTAGCATACTGTGGAGAACCATTA AAAATTTGAGAAGAAACAACAAGTATTATGTCACTTACTTCAAAGGCGTAGTTTTGGGAATTTGATGC AGTAAAGATTACCCTGTTTATGATTGTTCTTGAAGTCAAATGGGGACCTGTCCATTGTGCTCTAT TAATCTTGTGAGAAAAGTGTACCAAAAACAACTTGAAGTTGCTTGTCTAGGAGTTACTGGGTAG TTGTAAGTATTATTTTATTAATATAATGTAAAATAAAATGTTAAGTACTTAGTTTTGTTTTTCAAA GTAAAGCTGTAGTCAGCCTTATGTATGCCATTGACTCTGAAATGTATACCAGCCTTCACTGTGTACCG TGTGTATATAAATCCACAGAACCGGATGAGCTGCTTAGGGAGGGAATATTTCAAAGTGTACCAAGGAC CAAATCTGGAGTTCTCCCACTTTAGAGGATGGAAAGGGGCAGAGTAATCTAGCAAAGGAGACTGAGG CCAGTGAAGTAGGAAGAAAGTATTTCAAGGAGAGTGATGATTCTGTGAATATTGCTGAGAATTCAAATA AAAAGAGGACTGAGAAGTGA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001166017.2</u>
Summary:	This gene encodes a component of the spindle and kinetochore-associated protein complex that regulates microtubule attachment to the kinetochores during mitosis. The encoded protein localizes to the outer kinetochore and may be required for normal chromosome segregation and cell division. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Oct 2009]

Locus ID: 221150

MW: 57.9