

Product datasheet for **SC215057**

WSB2 (NM_018639) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	WSB2 (NM_018639) Human 3' UTR Clone
Symbol:	WSB2
Synonyms:	SBA2
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_018639
Insert Size:	1539 bp



[View online »](#)

Insert Sequence: >SC215057 3'UTR clone of NM_018639
The sequence shown below is from the reference sequence of NM_018639. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon **Red**=Cloning site

```

GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
AAAGAGTTCCTCACATACAGGACTTTTAAAGCAACACCACATCTTGCTTCTTTGTAGCAGGGTAAAT
CGTCCTGTCAAAGGGAGTTGCTGGAATAATGGGCCAAACATCTGGTCTTGCAATTGAAATAGCATTCTT
TGGGATTGTGAATAGAATGTAGCAAAACCAGATTCCAGTGTACTAGTCATGGATCTTTCTCTCCCTGGC
ATGTGAAAGTCAGTCTTAGAGGAAGAGATTCCACTTGCACGGCAACAGAGCCTTACGTAAATTTTCAG
TCCAGTTATGAACAGCAAGTGTTGAACCTCTTCTGCTTGTGTTTGAATCAAAGTGCAGTTACTGATGTTG
TTTTGATTATGCAACTAAGTAGGCTCCAGAGCCTCTCTAGTGGCAGAGCAGCTCACACTCCCTCCGCT
GGGAACGATGGCTTCTGCCTAGTACCTATCCTTGTGTTTCTGATGCAGTGGTAGCATTGGTTCAAGTTC
TCTCCTGCTGTGGTCAGAGTTGCTTCGATGTTGGCCAAGTGCTTTTCTTCTGGGCTCCCTTCTGACCT
GCAGGACAGTTTTCTGGAGCCATTTGGTATGAGGTATTAATTTAGCTTAACTAAATTACAGGGGACTC
AGAGGCCGTGCTCCTGACCGATCCAGACACTATTACTGGCTTTTTTTTTTTTTTTTAAACATGGTGTG
CATGTGCAGGAAATGACAAATTTGTATGTCAGATTATACAAGGATGTATTCTTAAACCGCATGACTATT
CAGATGGCTACTGAGTTATCAGTGGCCATTTATTAGCATCATATTTATTGTATTTCTCAACAGATGT
TAAGGTACAACGTGTTTTCTCGATTATCTAAAAACCATAGTACTTAAATTGAACAGTTGCAAAGATG
TCTTAATTGTGTAAGAATTGGTGTAGTCATGACTTTAGCTGATACTCTTATGTACGAGATCTGTCTCT
GCTGTTTAACTTCATTGGATTAATCAGCTGGTTCAACTCTACTGCGAAACAAAAATAGCTCCTTAAAA
GTACTGTTCTCCTTCAGTGGCATGTAGTTATCTAATCAAGACACCTCATTCAAACAAAACCTGCCTTAG
GAAAAATTAATATATTTTAAATTTTAAAGAAATACAACATCTTATTCTTTAGCTTTCTTAATCGG
TGCTTTATGGAGGCCAGTGAACGTTACATGACTCGTTGAGAAAGTTGAGGAATTCCTCTACCACCTT
TGTTGCTTGAAGAAAAACATGTCTTTTCAAAATGAGAGGCTTTCATTGAAGAAAAAGAAAAACAACAG
TTAAAAGCTTTTGGCTCTCTGTTTCATTTTTTCCATTAAGAAAAAAAAGTCCCTTTTAAAAACAAG
CAGGGCTGGGTGTGACGGCACATGCCTGTAAATCCCAGCACTTTGGGAGGCTGAGGTGGGAGGACTGCTT
GAGCCCAGGAGTTCAAGATCAGCCTGGGCAACAGAGTGAGATAAGACCAATCTTTATTTAAAAAATAA
ATAAATAATAGCCAAAAAAA
ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACGCCGCTTCTATGAAAGG

```

Restriction Sites: Sgfl-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: NM_018639.5

Summary: This gene encodes a member of the WD-protein subfamily. The encoded protein contains five WD-repeats spanning most of the protein and an SOCS box in the C-terminus. The SOCS box may act as a bridge between specific substrate-binding domains and E3 ubiquitin protein ligases. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jun 2013]

Locus ID: 55884

MW: 59.6