

Product datasheet for **SC215067**

FAM86B1 (NM_001083537) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	FAM86B1 (NM_001083537) Human 3' UTR Clone
Symbol:	FAM86B1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001083537
Insert Size:	1530 bp



[View online »](#)

Insert Sequence:	>SC215067 3'UTR clone of NM_001083537 The sequence shown below is from the reference sequence of NM_001083537. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAGCGATCGCC GAGATGGCAATGCTGAACCTCACACTGTAGGACTCACACATGACTCCAAAGGGATTGTGAGAATCAAGT CACTCTCGTGGGAAGAATTTTTATATGGGAAAGTGGATAAACTTTTCATTGGACTGGAATGTTTGGAGA TTATTAACTCCAAATCAGGAATCAAACTGCCCTCTAATAAGACATTAGCTGTCTAGGCGTGTGGGTG CCCCCTTTCTGCCAGCAGTTCTGGTTCTCAAGAAAATCATATAATCAGACATGAAAATCTAGCTCC AAAAATAGCATTTTCTTTCTGCAAATAAAAACGTGTGTATCAAGGATGACGTTCCCCCAACGTGGACAC ACTCGTTCTCTAGAAAGCCAAGCCCGTGCAGCTGCCACATCCCTGGACACACTCGTTCTCTACAAA GCCAAGCCCGTGCAGCTGCCACATCCCTGGACACACTCGTTCTCTACAAAGCCAAGCCCGTGCAGC TGCCACATCCCTGGACACACTCGTTCTCTACAAAGCCAAGCCCGTGCAGCTGCCACATCCAGGGCT TATGGTGCAGCAGCGCTTTTTTCAAGACAGGAATGAAAGTGTTAGGAACACGGCAGAAAGTGACACC TGGAGACCAATGCAGGATGAGGAGTACTGCAGAGGTACAGGGAAGTACAGAACAGTAATACGCTAG CAGGGGCATGGGCGATGAAGAACAGAAGACAGGAAGCATTTTCAGAGACTCCAAAGAAGAAATCAGGGCC AACCACAGCTTCCCGGTTCATTCACCAGGTGACACCACTGCCGTATTTAGCTTCTGGCCACTGGGAG GCGCTGCTCGAAAGGCTTTGCCCTGAGACTCCAAGAAGCTGCGGGAAGGACAGCAGGGGCCCTGGGGT TTAGCCTCTGGCCAGGAGTTATGTGTCCATAACCAAGGGAGCAGCTGCACCCAGCTCTATCCC ATCAGAGCTGCTGAGACTCTTGCAAGTTCTTCCAGAACTGGCTTAGCTTGCTGCAGGATCAGAAAAGT TTGAGAAAAGCATCTGCAAAATGCTAAAGAGCAGAGCTTAATCATTGCCTGTCCCCACCCATCCCAG GTCACCACTTGGATGACCCAGGTCCCCGACCCACCAACAACCCCTCCCAAGTCCCTAACTCCCTCACT TGGACTTGAGACCCTTCAACCCCCCTCACTTGGACTTGAGACCCTTCAACCCAGCAGCGCTCCAC TTCCAACCTGACATCATGCTTTCTGAAACTTCCCCGTATGTCCCACTTTCCCACTTGGTGCCCTGG AGCACCTTCCGGCCTCTACATGCTGTACGTTCCCTGTGAGCACCTCTCTCGGCCTCTGGCCAAACAC AGTCCCACCCATCTGTGGGTAACAAGGGGGTGTGGGTGTTCTTTTCAGCCTTGCTAAACTGTCTGAATC AAGGATCAGAACTACAGCCTGCAGGCCAAATCCAGCCACAGCCTGTGTTGTAATAAAGCTTTGTT GGAACAAAGCCA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001083537.4</u>
Locus ID:	85002
MW:	55.8