

Product datasheet for **SC215085**

ATP6V0A1 (NM_001130020) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ATP6V0A1 (NM_001130020) Human 3' UTR Clone
Vector:	pMirTarget (PS100062)
Symbol:	ATP6V0A1
Synonyms:	a1; ATP6N1; ATP6N1A; Stv1; Vph1; VPP1
ACCN:	NM_001130020
Insert Size:	1489 bp



[View online »](#)

Insert Sequence:	<p>>SC215085 3'UTR clone of NM_001130020 The sequence shown below is from the reference sequence of NM_001130020. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site</p> <pre>GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC CATATTCGGGAAGGGAAGTTTGAAGAGTAGTCCCTGTGAGGGCCGTGTGCCCATGCTACCTCCCCG CCTCCCTCCACAGTGATCAGCTGTGCCTCTCTGCCTGTTGGTTGTGATCTGTGGGACCAGCTCATTTCG TGTACCCTGTCTGTGAGTCATTTAGATAGAATAGTCCTCCTTGGGTCTCCACCACCCCTAGCTTTGT GTGTAGTGTAGTATTTTCTGGCTGTCACTCACTACTGAGCCAGCCTTGCCTCTTAGCCTCCA TCCATCCAGACAGCCCTCCACCTCCTGGTGGTGAAGCCAGTCTGCATCCACGCCATCCCAAAGCCC TTTCATCTCCCCGTGCATTGTAGATGGAAGGAGACCCATGCCATTACATCTAGACTTTGAGTCCC TGCATCTGCCACCGTAGTTTCTAGCAGGAGTAGTGGGGGAGTAATACAGATTCTCCCTAGAAGGGGA CACTGGTAACATGTCCACTCTTGGATTAGCAGGGGTGGGTCCAGGAAGATGATATTTGCGTCTTTTGC CCACCCCTGGCATTAGCTGGACCCAACTAGGCCATCATGAGTGGCTTCTCCCTGTCATCCCCAGGG GTCATAGGATATCTACACCGCCTTTCTGACCCACCCTGCACTCCCATCCTTTCTCTCTCCCCGTTC TGCCCTGCACTACATAGCACAGCCGGGATGCTTGGAAACAGAGGCCTTGGCTGCTCCGAGTGACAGGG CTTCCCTCTCTCGGGGTTGGCTTCTTCCAGGCCTTGCATGGGCCTGCCACAAGCACACCCTCAGGC CGAGGGTGCAGACTGATGCTCTTCCCTGATGGAGACCCTGAGATCTTCCCAACCCCAATCATGATGTC TTCAGTGTGGGACTGGGGTCTTGGTTCTGCCTGCAGCCTGCCTGGCTCCGCCCTAGTGCCTCCCTC CTCACCACTGGCCCCAGGTCTCAGGAGGGGTGCTCTGGGCAGGGAAGGTCAGTGTCACTGATGGTTT GCTGTTTGAAGCCATTGGCAGGGCTCCGTGCATGTGGCTGTGAGGGCTGCACAGTCCCTGCCAAGGGG CTTCTCCTTGTACCCCGAACCTTGAATCGTGTGCTGGCGTGGCAGCCCTGGCTAAGTTAATCCCA CCGCTTTCAGTGGTAGAAAGAATTCCCTGAGTGGCCAGGCTGGTGCCCTCCTCCTACCCTGGCTTTTC TGAGTGAGCTGCCTGGAGCCCTCATCCCTCTCCAGGCTGGGCTGGCCCTGGGCGGGGCCACTGTGTG CTGGCCACTGTGACCTGACCCGACCTTGTGCAGCCCCCTGCCCTGGTGTCTGGGTTTTCTGTGATGA TCTTTGCTCTGTTCCAGTGGGTTTTGAAGCAGAGTTCAGGGAACCCTGCCCAAGGTCTCCTGTTCAG ACATTCTATGTTGAATAAAGTATGTTTGACTTCCCCGGA ACGCGTAAGCGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCCAACCTGCCATCA CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG</pre>
Restriction Sites:	Sgfl-Mlul
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	NM_001130020.3

Summary:

This gene encodes a component of vacuolar ATPase (V-ATPase), a multisubunit enzyme that mediates acidification of eukaryotic intracellular organelles. V-ATPase dependent organelle acidification is necessary for such intracellular processes as protein sorting, zymogen activation, receptor-mediated endocytosis, and synaptic vesicle proton gradient generation. V-ATPase is composed of a cytosolic V1 domain and a transmembrane V0 domain. The V1 domain consists of three A and three B subunits, two G subunits plus the C, D, E, F, and H subunits. The V1 domain contains the ATP catalytic site. The V0 domain consists of five different subunits: a, c, c', c", and d. Additional isoforms of many of the V1 and V0 subunit proteins are encoded by multiple genes or alternatively spliced transcript variants. This gene encodes one of three A subunit proteins and the encoded protein is associated with clathrin-coated vesicles. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Locus ID:

535

MW:

54.2