

Product datasheet for **SC215094**

CERKL (NM_201548) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	CERKL (NM_201548) Human 3' UTR Clone
Symbol:	CERKL
Synonyms:	RP26
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_201548
Insert Size:	1542 bp



[View online »](#)

Insert Sequence: >SC215094 3'UTR clone of NM_201548
The sequence shown below is from the reference sequence of NM_201548. The complete sequence of this clone may contain minor differences, such as SNPs.
Blue=Stop Codon Red=Cloning site

```
GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG
TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC
GGAAGCATGGAAGAAATGATTCCAAAGTAAACATGTAATTGTTTCTAAAAGAAATGTACAACTAGAATT
TTAATATAAAGAAATGTGGCCATGCTTTAAACCTGTGGTCCCCACCCCTTTTGGGACCAGTTTCATGGA
AGACAATTTTTCCGTGGATGGGGGAATGGTTTCAGGATGAACTGTTCCACCTCAGATCATCAGGCATTA
GATTCTCTTGAGGAACATGGCACTCTAGATCCCTCACATACCAGTTCACCATATGGTTCGCACCTCTAT
GAGAATCTAATGCTGCTGCTGATCTGACAGGGGGAGGAACTCAGTAATGCTCCATTGCTCTAACCTCCT
GCTGTGTGGCCAAACGGATATCCTAATAGAAATGATGTTACTAACAAATAGCCTCTAGGCTTTAAAA
GTCAAGATGAAATATCAATTCAAATCTGTTTTAGATATCAATTCAAATCTTTTTTTTACAATGTGT
TTATTGGACACACAAAAAACTTTGCAACCATATAATACATCAATTTAACCTAGATAATTCTGAAA
TAATTTGGATTCTTTCATTTTTTCAGGATTTGAGCTCATCAATTATGTTAAATTTCTATATTCTGTTAC
AAATATAATACAGATTTTATAAGTCTGCCTTGATTCAGTCTCCCTAATCCAGAGCAGTTAAATTTCAA
AATCACATAGTTATATAGCCTACCTGCCAAAAAGAATAACCCATCCCTGTCTTCTTAAAGTGGTAGGG
TTTTAAAATTCATGATACAAGTAACAAAGTAAGTTGATAACAAACTACTATTTTGAAGTGGCAAAATCT
AGCCCAATATGTAATTGTGTTCCCTTGTATCTGAGAGGCCTGTGTTGGACAACCATAGTACTCTGCTT
TCACTGTGCTCAGATAACCTGAAATATATAGTTGTTTCTGAGTAGTCTGAGTAGGAAAGAACCAATG
TAAATTTCTTAGAGACCATCTGGAGAAGACAGTCTAAATCGGGGCTTGATTTTTTTCATCTATGATCC
TAGATGTTCTCTGAAAGAATGACTAGCCTTTTTATAACGATACAAAGTTTTTTAAAAATACATAAATG
TCATAAAGTTTTATCCTGTGTATACACACTCTGGGTTCCAGGAGTTTTGTGCTTACATTTCACTACTCA
AATGCCTTCGCAATCTCAGAACATTCCTCACCATCACTGCAGGCCTGTGAACAGATAACTTTCATCAA
TTAGGATGGCTGGCTAGTGGCAGGGACTTCATATAAAGAAATCTTTCCAGAATTGTTTTCCACATTG
CCTTCTAATTTGATCACTCCCTATAACAGATGTTCAAAATTCGGCCACAGCTCTTATGAAGGCAGAG
TCAGATAAAATGATATTCTGCAATTATCTAGATCATAAGCATAAAATAAATTACAAAAATGATTGTA
TTTTAAAGAAATTTTTAAAAATCCAGAAAGTCATTTAAAAATAGAACCTCATATAGTATGAACAACTATA
AAAAATTTTTACATTTTCATGAAAA
ACGCGTAAAGCGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA
CGAGATTCGATTCCACCGCCGCTTCTATGAAAGG
```

Restriction Sites: SgfI-MluI

OTI Disclaimer: Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences, e.g., single nucleotide polymorphisms (SNPs).

Components: The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.

RefSeq: [NM_201548.5](#)

Summary:

This gene was initially identified as a locus (RP26) associated with an autosomal recessive form of retinitis pigmentosa (arRP) disease. This gene encodes a protein with ceramide kinase-like domains, however, the protein does not phosphorylate ceramide and its target substrate is currently unknown. This protein may be a negative regulator of apoptosis in photoreceptor cells. Mutations in this gene cause a form of retinitis pigmentosa characterized by autosomal recessive cone and rod dystrophy (arCRD). Alternative splicing of this gene results in multiple transcript variants encoding different isoforms and non-coding transcripts. [provided by RefSeq, May 2010]

Locus ID:

375298

MW:

59.4