

Product datasheet for **SC215104**

E74 like factor 1 (ELF1) (NM_001145353) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	E74 like factor 1 (ELF1) (NM_001145353) Human 3' UTR Clone
Symbol:	E74 like factor 1
Synonyms:	EFTUD1; RIA1
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001145353
Insert Size:	1557 bp



[View online »](#)

Insert Sequence:	>SC215104 3'UTR clone of NM_001145353 The sequence shown below is from the reference sequence of NM_001145353. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site <pre> GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAGGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AACGAACTGCTGGAACCAACTCTTTTAGTTAATATACCAAAGCTTATGAATAATTGTTTGTAAATTG AACATTTTCAATTATATGCAGACTGACTGATTCTAAGATAAATTCTAAGGAGTTTCTAATTTTGTAAAT TGTTAAAAATAGAGTTAATTTTGACTTTGTAGATGAGGGAGGAAAACCTCAACTGTTTCTCTTTGTAT CTAAATGTTTCAGAATTCATCGTGAAGGAACAGGCATTTTACACTATGAAGACATTCTTTTGAGATTT TTATTTTCAGTTGCTATATCATAAGCATTTTTAAAGTTTCTTTTCTAATTTTACATTGTATTAGATTTTC TGATTCTTTTGTAAATACAGAACTTAAATAGAAGGCAACAGGAAATTTATATAGGAACTATTTTCATTC CACTTGCTGAAGTTAAGTCTTGACTCTTCAAATGCAAAAAACCTATTTTATGCTTTGTTAAATATG GTGTCACCTTAGATTGACTTTAGTTGACTGCACTATATAATATAGAAGTATGAATATGTAGAATAACATG AAAAATTGGAGGTGCTGGTGGTATGGCTGACCCTGTTTCAGAAGCAGGATAGTATAAAGCATCAGCCT AAGAATGGCACTCCCACTAACTAGCTATGAATCTTGACCTCTTTGGGCTTTAGTTCTCTCATAAAAG GAAGAGATGTATTGGATTAGACTAGATGATCACCCTTTCTCTTAGTTCTAATTTTTTAAATTCTAA TACCTATATTTTCAAGTTATGTCAATTAATCATTATCAGGTTATTTCTAATGTAAAGATAGCTAAAA TGTTGCAGAGAAATAAGTGACCAACAAAAATTTATTCATCTGTTATGGGTAAGATCTGCCATAAATTCT TCCTAAATAATTTGTTTACTAACTCTTTAGGCCACTGTGCTTTGCGGTCCATTAGTAACTTGTTGTC TAAGTGCTAAACAGAATACTGCTATTTTGAGAGAGTCAAGACTCTTTCTTAAGGGCCAAGAAAGCAACT TGAGCCTTGGGCTAATCTGGCTGAGTAGTCAGTTATAAAGCATAATTGCTTTATATTTTGGATCATTT TTTACTGGGGGCGGACTTGGGGGGGTTGCATACAAAGATAACATATATATCCAACCTTTCTGAAATGAA ATGTTTTTATAGTTACTTTTTCACTGTAATAATGTACATTTAATGTACAAAGAAAAAATGTCTCTG CAAATTTTCTAGTATAACAGAAATTTTGTAGATGAAAAAATCATTATGTTTAGAGGTCTAATGTAT GTTTTCATATTACAGAGTGAATTTGTATTTAAACAAAAATTTAAATTTTGAATCCTCTAAACATTTT GTATCTTTAATTGGTTTATTATTTAAATAAATCATATAAAATTTCTCAGTGTCTGTTTTCAGGCAAAAGT TTCTTAAAGAATAAGTGTGCAGAGAATATTACTAGAACATCAGCATTACTTAATGTTTATAAATAAAT TCATTAGTCAGAATTGCAAAAAAAAAAAAAAAAAAAAA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG </pre>
Restriction Sites:	Sgfl-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001145353.1</u>
Summary:	This gene encodes an E26 transformation-specific related transcription factor. The encoded protein is primarily expressed in lymphoid cells and acts as both an enhancer and a repressor to regulate transcription of various genes. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Feb 2009]
Locus ID:	1997

MW: 61.3