

Product datasheet for **SC215114**

ZNF761 (NM_001008401) Human 3' UTR Clone

Product data:

Product Type:	3' UTR Clones
Product Name:	ZNF761 (NM_001008401) Human 3' UTR Clone
Symbol:	ZNF761
Synonyms:	ZNF468
Mammalian Cell Selection:	Neomycin
Vector:	pMirTarget (PS100062)
ACCN:	NM_001008401
Insert Size:	1543 bp



[View online »](#)

Insert Sequence:	>SC215114 3'UTR clone of NM_001008401 The sequence shown below is from the reference sequence of NM_001008401. The complete sequence of this clone may contain minor differences, such as SNPs. Blue=Stop Codon Red=Cloning site GGCAAGTTGGACGCCCGCAAGATCCGCGAGATTCTCATTAAAGCCAAGAAGGGCGGAAAGATCGCCGTG TAACAATTGGCAGAGCTCAGAATTCAAGCGATCGCC AGACTTCATACTGGAGAAAAACAAGTGTAAATGAATGTGGTGAGGTTTTTAATCAACAAGCACACCTTGC AGGTCATCATAGAAATTCATACTGGGAGAAACCTTAGAAATGTGAAGCATGTGATAAAGTTTACAGTGG CAAATCAAGCCTCAGAAGACAGGAGAATTCATACTGGAGAGAAAGCTTATAAATGTGAAGAATGTCACA AAGTTTACAGTCGCACATCAAAACGTGAAAGACAGGAGAATTCATACTGGAGAGAAACCATAAAAATGT AAGAGTTTGTGACAAGGCTTTTGGGCATGATTCGCACCTGGCACAACATGCTAGAATTCACACTGGAGA GAAACCTTACCAGTGTATGAGTGTGGCAAAGCCTTTAGTAGGCAGTCAACACTTGTTTACCGTCAGGC AATCCATGGTGTAGGGAACTTTACTAAGGTAATGATTGTCACAAAGTCTTCAGTAATGCTACAACCAT TGTGAATCACTGGAGAATCCATAAGGAAGAGAGATCATACTAGGGTAATAAATGTGGCAGATTTTTTCAG ACATTGTTTCATACCTTGCAGTTCATCGGTGAACCTCAACGCTGGAGAGAAACCTTACAAATGTCATGACT GTGGCAAGGTCTTCAGTCAAGCTTCATCTATGCAAAACATAGGAGAATTCATACAGGAGAGAAACCTC ACATGTGTGATGATAGTGGCAAAGCCTTCACTTCACACCTCATGAGACATCAGAGAATGCATACTGGAC AGAAATCTTACAAATGTCATCAATGTGCCAAGGTCTTCAGTCTGAGTTCCTCCTTGAGAATATGAGA AAATTCATTTTGGAGGTAGTTGGTCCATATGCAATGAGTAGAGCAAACCATCAAGCATTAAATTGACATT AGGGTCAATTCAGCATTGACTTGAGTTTGTATTGACTTAACATTGAGTTCAGCATTAAATTGACATTAG TGTTTATGTAAAGAGGATTGGGCCAGGCACATCAGCTTACACCTGTAATCTGAGCGCTTTGGGAGGCCA AGGTGGGTAGATCACTTGAGGTGAGGAGTTTGAATCAGCCTGGCCAACAGACGTGAGCCATTTTCCCA GCCTGTTTTTTTGTCTTTAAAAAACTGATAGGGATTTTATGGATATCATGTTGAATCTAAATCACA TTGGGTTATTATATAATCATTTTCAATATTAATTTTTTCAAGCTATCAATATGGTTGTAGCTCAATG TTTTTAATCATTTTGATCAATGTTTGTAGATTTCAAGGTACAACTTCTGACCTTTGTACGTTTATTTTC TAAGTATTTCTTAAGTTCTCCAGCAAATGGAAGTGTAAAAATTTCTTTTAAATTTGTTTATTGTT AAAGTATGGAAATTAACATAATTTTGGTGCTGATACTGTATTGTGCAATCCACTGACTATGTTACTT AGTTCCAGTAGTATTTTGGTTGGCTCTTTGTGATTTTCTACACAGAAGATTATGTCATCTACAAACAAA TATAATTTTACTTCTTACTTTCTGA ACGCGTAAGCGGCCGCGCATCTAGATTGAAGAAAATGACCGACCAAGCGACGCCAACCTGCCATCA CGAGATTTTCGATTCCACCGCCGCTTCTATGAAAGG
Restriction Sites:	SgfI-MluI
OTI Disclaimer:	Our molecular clone sequence data has been matched to the sequence identifier above as a point of reference. Note that the complete sequence of this clone is largely the same as the reference sequence but may contain minor differences , e.g., single nucleotide polymorphisms (SNPs).
Components:	The cDNA clone is shipped in a 2-D bar-coded Matrix tube as 10 ug dried plasmid DNA. The package also includes 100 pmols of both the corresponding 5' and 3' vector primers in separate vials.
RefSeq:	<u>NM_001008401.4</u>
Summary:	May be involved in transcriptional regulation.[UniProtKB/Swiss-Prot Function]
Locus ID:	388561
MW:	59.2